



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## X-Ray/Radiographic Positioning Doll

[Read More](#)

**SKU:**

**Categories:** Clinical Skills Trainers



### Product Description

This life-size Radiographic Positioning Doll provides a safe, realistic training tool for radiography students. With a fully articulated, radiolucent skeleton and identifiable anatomical landmarks, it enables accurate practice in patient handling, positioning, X-ray techniques, and CT scanning. Durable, lightweight, and supplied with individual X-rays for reassurance and authenticity.

Durable and lightweight, the AR10A is designed for repeated use in teaching environments. Its realistic design helps build assurance in technique and accuracy, while also encouraging correct use of equipment and positioning methods. By bridging the gap between classroom learning and patient contact, the Doll allows students to refine their skills in a safe, controlled setting, reducing anxiety and improving readiness for clinical practice. The opportunity to make and correct mistakes without consequence greatly enhances understanding and retention, making it an invaluable step between theory and patient care.

Adopted by institutions worldwide, it continues to be a trusted training resource for developing skilled and capable radiographers. Its proven durability and adaptability ensure long-lasting value, supporting both individual learning and group demonstrations across educational settings. Institutions consistently report improved engagement, as the realistic structure encourages hands-on participation and reinforces good practice from the very beginning of training.

## Skills

- An invaluable teaching aid in Schools of Radiography, without bringing discomfort or danger to patients or students
- Widely used by Schools of Radiography throughout the world
- Life size human figure for training Radiographers in the following skills:
  - Positioning techniques
  - Patient handling
  - X-Ray machine use

## Features

- Life-size figure containing a fully flexible articulated skeleton produced from a specially selected grade of plastic
  - No metal parts are used in the articulation of the skeleton
  - Anatomical landmarks easily identified
  - The body covering has the combined properties of being able to simulate surface anatomy whilst remaining totally radiolucent
  - Representations of the larynx, heart, lungs and kidneys are constructed from radiolucent material
  - Weight – the doll is light enough to be carried easily
  - Joint flexibility – there is slight hyperflexion of the knees and elbows.
  - The hip rotation is sufficient and generally the joints are realistic
  - Flexible and durable articulation of the skeleton ensures a lengthy period of use
- 
- Height: 165cm

- Weight: 12kg

## **2 Year Guarantee**

- This model, manufactured by Adam,Rouilly, comes with a 2 Year Guarantee. This guarantee applies to models which have been used correctly and covers durability and functionality.

## **Individual X-Rays**

### **Individual X-Rays for each Doll**

- Each Doll is individually X-rayed
- Customers can be reassured that there are no metal parts in the Doll by referring to the set of digital X-Rays which are supplied with each Doll
- Images shown were taken using medium frequency generator; no grid; fine focus; 100cm S.I.D. Fuji CR system

## **Positioning, CT & MRI Scanning**

### **Positioning**

- The Doll is best used in the recumbent positions as, although it is possible to use an upright bucky with supports, this makes the use of positioning blocks and immobilization devices more difficult
- The Doll will lie naturally in the neutral position, prone or supine, without support
- A wide selection of positioning aids must be at hand, e.g. foam blocks, wedges, sand bags and compression band. These are essential to maintain any position other than neutral as the skin of the Doll has some resistance to overcome. This is, of course, an advantage in the teaching as it prevents the student cutting corners by omitting the use of these positioning aids
- The Doll can be positioned for all standard projections with no more difficulty than a difficult patient might present

### **CT & MRI Scanning**

- The X-Ray Positioning Doll contains no metal parts so can be imaged for positioning and centering for whole body CT and MRI Scanning
- Good image quality in terms of body outline and fundamental structure can be obtained. Internal

organs are less clearly defined.

- Certain areas of the skeleton structure e.g. thoracic vertebrae and the head of the femur gave similar values to the typical patient
- As the Doll is very different from the human body in terms of attenuation it cannot be used for training in, or practising exposure control

## **Manufacturer**

Adam,Rouilly