



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Intubation Head for CPR

Item No. 1019711 [P71/AH]

Weight 1.6kg

[Read More](#)

SKU:

Categories: Airway Management Adult



Product Description

Upgrade CPRLillyPRO™ to include airway management training with the Trucorp® advanced airway head: The head offers the most realistic internal anatomical details for a full range of airway management techniques. It is easy to install on your existing CPRLillyPRO training manikin and can quickly be interchanged when switching between CPR and airway management training scenarios.
Life-like airway for realistic skill training in:

- Double naso-tracheal intubation
- Bag and Mask ventilation techniques
- Full range of supraglottic devices
- Direct laryngoscopy
- Endotracheal tube insertion
- Awake fibre optic examination
- Combi tube insertion

Recommended equipment size for airway management training:

- 7.0-7.5 mm ID for nasal intubation
- 8.0-9.0 mm ID for oral intubation
- Size 3-5 for LMA laryngeal masks
- Similar respective sizes for other supraglottic devices

Manikin can be used with Supraglottic airway devices for airway training. Allergy note: head skin contains latex.

Why Choose the TruCorp Intubation Head for CPRLillyPRO™ for Your Simulation Lab in India?

The TruCorp Intubation Head for CPRLillyPRO™ (Item No. 1019711, MPN: P71/AH) by TruCorp®, Ireland is a 1.6 kg advanced airway training head designed as a direct upgrade for the CPRLillyPRO™ training manikin — allowing simulation labs that already own the CPRLillyPRO™ to add the full range of advanced airway management training capabilities to their existing CPR platform without replacing the manikin. The head installs quickly and can be interchanged between CPR and airway management configurations, enabling instructors to transition the same manikin between CPR-only sessions and integrated CPR-plus-airway training sessions without disruption. For simulation labs in India that have standardised on the CPRLillyPRO™ for BLS training and want to extend into ACLS airway management without purchasing a separate dedicated airway trainer, this is the correct and most cost-effective upgrade path.

The TruCorp advanced airway technology delivers the most realistic internal anatomical details available for the full range of advanced airway management techniques — the same airway quality that defines the TruCorp AirSim® range. The procedure set covers: double nasotracheal intubation, BVM ventilation techniques, the full range of supraglottic device insertion, direct laryngoscopy, ET tube insertion, awake fibreoptic examination, and Combitube® insertion. Recommended device sizes are 7.0–7.5 mm ID for nasal intubation, 8.0–9.0 mm ID for oral intubation, and size 3–5 LMA laryngeal masks. The manikin is compatible with the full range of supraglottic airway devices for extended airway training.

?? Allergy notice: The head skin contains latex. Users with latex sensitivity or allergy should take

appropriate precautions before use.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of TruCorp® products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: Which manikin is this intubation head compatible with? A: This head (P71/AH) is designed as an upgrade for the CPRLillyPRO™ training manikin by TruCorp. It is not a standalone product and requires the CPRLillyPRO™ manikin for use. Contact us to confirm compatibility with your specific configuration.

Q: Does the head skin contain latex? A: Yes. The head skin contains latex. Users with latex sensitivity should take appropriate precautions before use.

Q: What is the warranty and how do I get support? A: TruCorp products carry a manufacturer warranty. For support in India, contact SEM Trainers at sem@semtrainers.com or +91-88495 63724.