



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Injection Belly

Item No. WA28395U
Weight 0.91kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

This soft stomach replica will be an invaluable tool for teaching patients how to properly insert and rotate the infusion set for insulin pumps, avoiding the 2" (5 cm) area around the navel, as well as giving self-injections of many types including pen-style injections. The trainer has lifelike skin with what feels

like real stomach tissue underneath. There is a soft plastic backing to keep needles from going through and holes to allow the material to breathe and dry, should you choose to inject liquids. Distilled water may be injected, but the more liquid that is injected into the simulator, the longer the time for evaporation. Actual insulin or other medications should not be used, as they may deteriorate the material. The injection belly is considered disposable; however, with proper care and treatment it will serve its purpose for a long period of time. Using the smallest needles possible when simulating injections and rotating the insertion points will help prolong the life of the trainer. Dimensions: 10" L x 6" W x 1-1/4" H.

Why Choose the Injection Belly for Diabetic Patient Education in India?

The Injection Belly (Item No. 1013056 [W44765], 0.91 kg, MPN: WA28395U) is a soft, realistic stomach replica designed for teaching patients how to properly insert and rotate insulin pump infusion sets — avoiding the sensitive 2" (5 cm) area around the navel — as well as practicing pen-style and other self-injection types.

The trainer has lifelike skin with realistic underlying tissue feel, backed by a soft plastic layer that prevents needles from passing through, with holes that let the material breathe and dry after liquid injection. Distilled water can be used for practice injections, though larger volumes take longer to evaporate.

Important care notes:

- Use distilled water only — never actual insulin or medications, which will deteriorate the material
- Use the smallest needles possible and rotate injection sites to extend the trainer's life
- Considered a disposable product, but will serve its purpose for a long period with proper care

Dimensions: 10" L × 6" W × 1¼" H.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of leading diabetic patient-education equipment in India, with over 30 years of experience supplying medical training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Can insulin be used for practice injections? A: No — only distilled water should be injected. Actual insulin or medications will deteriorate the material.

Q: Is this a reusable product? A: It's considered disposable, though with proper care (small needles, rotating sites) it serves its purpose for a long time.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.