



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Willi and Cathi Catheterization Trainers

Item No. SB50115
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:Catheterization



Product Description

Trainers fit into Interchangeable Catheterization & Enema Simulator, (LG02099U">LG02099U, In-Service Home Care Training Simulator, (LG02098U">LG02098U, and most Laerdal® manikins or into the remedy stand (coming soon). Material feels realistic and delicate like actual tissue. Use the Cathi female catheterization trainer to practice locating the urethral meatus and for urinary catheter placement. The labia minora realistically closes over the vaginal and urethral openings. The tissue is supple and pliable for easy spreading with thumb and index fingers. The Willi male catheterization trainer has a removable foreskin and can simulate an enlarged prostate. Successful catheterization on either trainer will be confirmed with urine flashback. Manikin not included.

Full page content retrieved. Lean OEM description — one paragraph only — but enough to build from. Here is the complete audit and ready-to-paste copy:

Category: ? Correctly filed — Medical Simulators ? Clinical Skills Trainers ? Catheterization. Correct.

Item No. **SB50115**, Make: **Nasco Healthcare, USA ?** — Make/Brand field already populated. **Weight field is blank** — needs to be obtained from Nasco Healthcare datasheet. Description is a single OEM paragraph with no India-market framing, no structured benefits narrative, no clinical context, no FAQ, and no meta description.

Important structural note: This is a **two-in-one listing** — Willi (male) and Cathi (female) are two separate trainer inserts sold together under one SKU. The description must make this clear. The urine flashback confirmation feature is the key functional differentiator. Compatibility with three manikin families (LG02099U, LG02098U, and most Laerdal® manikins) is a critical purchasing decision point that the current description underplays.

Note on the title: "Willi and Cathi" are Nasco Healthcare's registered product names for these trainers — retain them in the title exactly as written.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Willi and Cathi Catheterization Trainers for Your Clinical Skills Lab in India?

The Willi and Cathi Catheterization Trainers (Item No. SB50115) is a dual-gender urinary catheterization training set from Nasco Healthcare, USA — comprising two interchangeable anatomical inserts: **Cathi**, the female catheterization trainer, and **Willi**, the male catheterization trainer. Sold together as a single kit, they allow nursing colleges, medical schools, urology departments, and clinical skills labs in India to train both female and male urinary catheterization technique on the same compatible manikin platform — switching between genders in seconds for skills station rotation, OSCE setups, and comparative training sessions.

Both trainers are designed to slot directly into three compatible Nasco Healthcare manikin systems and most Laerdal® manikins:

- Interchangeable Catheterization & Enema Simulator — **LG02099U**
- In-Service Home Care Training Simulator — **LG02098U**
- Most Laerdal® adult manikins

Manikin is not included. The SB50115 kit contains the Willi and Cathi inserts only.

Cathi — Female Catheterization Trainer

The Cathi trainer trains the complete female catheterisation workflow: identification of the urethral meatus, correct labia separation technique, and catheter placement. The labia minora realistically close over the vaginal and urethral openings as they do in a real patient, requiring the trainee to actively locate and separate the labia with thumb and index finger before catheter insertion — the step most commonly performed incorrectly by novice practitioners. The tissue material is supple and pliable, replicating the feel of actual tissue, and successful catheter placement is confirmed by urine flashback.

Willi — Male Catheterization Trainer

The Willi trainer trains male urethral catheterisation including foreskin retraction and management of a simulated enlarged prostate — two anatomical variables that significantly affect difficulty and technique in a clinical setting. The foreskin is removable, allowing training both with and without a foreskin present. The enlarged prostate simulation prepares trainees for the resistance and trajectory change encountered when catheterising a patient with benign prostatic hyperplasia — a clinical presentation frequently encountered in male catheterisation practice. Successful catheter placement is confirmed by urine flashback on both models.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states in India. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Frequently Asked Questions

Q: What is included in the SB50115 kit? A: The kit includes two anatomical trainer inserts — Cathi (female) and Willi (male). No manikin, catheter, lubricant, or urine simulation fluid is included.

Q: Which manikins are compatible with these trainers? A: The Willi and Cathi inserts are compatible with the Nasco Healthcare Interchangeable Catheterization & Enema Simulator (LG02099U), the In-Service Home Care Training Simulator (LG02098U), and most Laerdal® adult manikins. Confirm compatibility with your specific manikin model before purchase.

Q: How is successful catheter placement confirmed? A: Both the Willi and Cathi trainers confirm correct catheter placement by urine flashback — exactly as in clinical practice.

Q: Can the Willi trainer simulate an enlarged prostate? A: Yes. The Willi male trainer includes an enlarged prostate simulation, preparing trainees for the resistance and directional change encountered when catheterising a patient with BPH.

Q: Is the foreskin on the Willi trainer removable? A: Yes. The foreskin is removable, allowing training both with and without a foreskin to replicate the range of patient presentations encountered in clinical practice.

Q: Can these trainers be used independently without a manikin? A: The trainers are designed to be mounted in a compatible manikin or simulator. A remedy stand is listed as coming soon from Nasco Healthcare for standalone use — contact SEM Trainers for availability.