



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Jaw Thrust Kyle 3 Year Old CPR Manikin

Item No. 100-2960
Weight 8 lbs.
Brand Nasco Healthcare

[Read More](#)

SKU:100-2960

Categories:BLS Child



Product Description

This jaw thrust capable pediatric manikin offers you another option for teaching health care professionals. You will find all of the same features of our regular Kyle manikin with the addition of this simple rugged construction technique for instructing the procedure of moving the jaw and tongue out of the way without moving the head. Two-rescuer technique allows you to ventilate the patient without lifting the jaw, although that maneuver also functions with JT Kyle. Manikin comes with three face pieces, three airways, and carry bag. Size: 34" x 11" x 7". Sh. wt. 8 lbs.

Why Choose the Jaw Thrust Kyle™ 3-Year-Old CPR Manikin for Your Simulation Lab in India?

The Jaw Thrust Kyle™ 3-Year-Old CPR Manikin (Item No. 100-2960) by Nasco Healthcare, USA is approximately 3.63 kg, 86 × 28 × 18 cm — a variant of the Kyle™ CPR manikin family that adds a specifically engineered jaw thrust mechanism, designed for healthcare professionals who need to train the jaw thrust airway opening technique on a 3-year-old child patient. The jaw thrust — displacing the

mandible anteriorly using two hands positioned at the angles of the jaw, without tilting or extending the head — is the correct technique when cervical spine injury is suspected and head tilt/chin lift is contraindicated. This is the only configuration in the Kyle™ family that supports this technique. For PALS training centres, paediatric emergency medicine departments, PICU in-service programmes, and nursing institutes in India where trainees must be competent in both head tilt/chin lift and jaw thrust airway management on paediatric patients, this manikin provides both techniques on a single 3-year-old platform. Supplied with three face pieces, three airways, and carry bag.

The jaw thrust mechanism uses simple, rugged construction that allows the jaw and tongue to be displaced anteriorly without moving the head — practising the technique that protects the cervical spine during airway management in trauma scenarios. The two-rescuer technique is explicitly supported: one rescuer can hold the jaw thrust while a second ventilates the patient without lifting the jaw, or the jaw-lifting ventilation manoeuvre can also be performed. All standard Kyle™ CPR features are retained: correct 3-year-old proportions for paediatric CPR compression hand placement, disposable airways with individual face pieces for hygienic multi-trainee use, and realistic anatomical landmarks for paediatric resuscitation training. The face pieces and airways are individually replaceable to maintain hygienic practice standards.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: What does the jaw thrust mechanism add over the standard Kyle™ 3-year-old manikin? A: The jaw thrust mechanism enables training in anterior mandibular displacement without head tilt or extension — the correct airway opening technique when cervical spine injury is suspected. Standard Kyle™ supports head tilt/chin lift only. Jaw Thrust Kyle™ supports both techniques on the same manikin.

Q: Does this manikin support two-rescuer CPR training? A: Yes. The two-rescuer technique is explicitly supported: one rescuer holds the jaw thrust while a second performs ventilation without lifting the jaw — a key training scenario for paediatric trauma resuscitation.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.