



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

SimScope® Tablet with Software

Item No. 1021559
Weight 4kg
Brand Cardionics
MPN: 718-3440

[Read More](#)

SKU:

Categories:Auscultation



Product Description

You can program the SimScope to play desired sounds from our vast sound library which is already embedded in the SimScope tablet. The tablet is part of the SimShirt System (1021552).

SimScope® Tablet with Software — Programmable Sound Console for the SimShirt Auscultation Simulator

The SimScope® Tablet with Software (Item No. 1021559, MPN: 718-3440, 4 kg) by Cardionics, USA is the programmable control tablet for the SimShirt Auscultation Simulator System (Item No. 1021552). It comes pre-loaded with Cardionics' full auscultation sound library, allowing instructors to select and play specific heart, lung, bowel, or bruit sounds directly through the connected SimShirt for real-time student auscultation practice.

Built specifically for SimShirt, not a standalone accessory. Unlike SimScope stethoscopes and patch kits that work across multiple SAM-family trainers, this tablet is a dedicated component of the SimShirt System — it's the interface instructors use to drive which sound plays and when, making it central to running a live SimShirt training session rather than an optional add-on.

For simulation labs running wearable, standardised-patient-style auscultation training (where the SimShirt is worn over a live or task-trainer setup), this tablet lets faculty switch conditions on demand during a session — useful for OSCE-style exams or live demonstration teaching where students need to identify a sequence of different findings.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of Cardionics products in India with over 30 years of experience supplying medical simulators and auscultation training accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: Can this tablet be used with SAM 3G, SAM II, or PAT manikins? A: No — it's a dedicated component of the SimShirt Auscultation Simulator System (Item No. 1021552), not a general-purpose SimScope accessory. Confirm your system setup with our team before ordering.

Q: What can instructors do with the tablet? A: Select and play specific sounds from the embedded sound library, letting instructors control exactly which heart, lung, bowel, or bruit finding a student hears during a live session.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.