



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Arrhythmia Television Interface

Item No. 1018150 [W49429]
Dimensions 0.001 x 0.001 x 0.001 cm
MPN: ADI-10

[Read More](#)

SKU:

Categories: Auscultation



Product Description

In classroom situations, where no ECG-monitor is available, this Arrhythmia Television Interface (ATI) allows connection of the Arrhythmia Tutor (W49007) to a conventional television monitor. The selected arrhythmia rhythms will be displayed on the monitor, additionally it is also possible to freeze the action or to add ECG grid lines. The unit is supplied with power cords, NTSC television connector cable, video connector cable, and a three-lead cable. The unit is suitable for 115 Volt. (For countries with 230 Volt, a transformer unit, which can be locally purchased, is necessary). Please be sure your television is capable to receive the NTSC signal. PAL televisions can only be adapted if they contain a video inlet

(CINCH) in which the video connector cable can be plugged.

Arrhythmia Television Interface — Display Rhythms on a Classroom TV

The Arrhythmia Television Interface (ATI) (Item No. 1018150 / W49429, MPN: ADI-10) connects the Arrhythmia Tutor (W49007) to a conventional television monitor for classroom teaching where a dedicated ECG monitor isn't available. Selected arrhythmia rhythms display on the connected TV, with the ability to freeze the action or overlay ECG grid lines for detailed review.

The unit ships with power cords, an NTSC television connector cable, a video connector cable, and a three-lead cable, and is built for 115V power.

Important compatibility note for India: This unit outputs an NTSC signal at 115V. For use in India, you will need a separately-sourced 230V step-down transformer, and your television must have a CINCH/RCA-style video input to accept the video connector cable — this is uncommon on TVs manufactured in the last several years. Please confirm both requirements with SEM Trainers before ordering.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and auscultation training accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: Will this work with a standard Indian television? A: Only if the TV has a CINCH/RCA video input and you use a locally-sourced 230V step-down transformer. Most modern TVs (HDMI-only) are not compatible — confirm with SEM Trainers before ordering.

Q: What does this connect to? A: The Arrhythmia Tutor (W49007).

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 to confirm compatibility and pricing.