


**SEM****SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724**Mobile 1:** +91 98791 03905**Email:** sem@semtrainers.com**Website:** www.semtrainers.com

Front Teeth, set of 3

Item No. 1005397 [W30510]
Weight 0.05kg
Dimensions 1.5 x 3 x 3.5 cm
MPN: Front Teeth for M167

[Read More](#)**SKU:****Categories:** Airway Management Adult

Product Description

for item number W30508 NOTE: These are replacement teeth for medical training simulators and are not for human dental use.

Why Choose the Front Teeth, Set of 3 for Your Simulation Lab in India?

The Front Teeth, Set of 3 (Item No. 1005397 / W30510) is the genuine OEM replacement front teeth consumable for the Endotracheal Intubation Simulator (W30508) by Nasco Healthcare, USA, weighing 0.05 kg at $1.5 \times 3 \times 3.5$ cm per set. The front teeth on the Endotracheal Intubation Simulator are a breakout consumable — they are designed to disengage when a trainee applies incorrect laryngoscopy technique, specifically levering the blade against the upper teeth rather than using correct straight-lift technique. This built-in consequence mechanism is one of the most effective feedback tools for

correcting the most common laryngoscopy error: the upper-teeth lever technique that causes dental damage in real patients. When the teeth break out, training has served its purpose — but only if the simulator is restored to full training condition with genuine OEM replacement teeth before the next session. The 3-pack format provides a practical stock quantity for simulation labs running regular intubation training sessions.

Important notice: These replacement teeth are designed exclusively for use with the Endotracheal Intubation Simulator (W30508). They are not for human dental use and must not be used outside their intended purpose as a medical training simulator consumable. Only genuine OEM replacement teeth are guaranteed to fit the W30508 simulator correctly and maintain the correct breakout force that makes the consequence feedback clinically meaningful. Using non-OEM substitutes risks incorrect breakout resistance — either too low (triggering without genuine technique error) or too high (failing to break out when technique is incorrect) — both of which degrade the training value of the simulator's feedback mechanism.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which simulator are these replacement teeth compatible with? A: These teeth are designed exclusively for the Endotracheal Intubation Simulator (Item No. W30508) by Nasco Healthcare. They are not compatible with other simulator models.

Q: Can these teeth be used for human dental purposes? A: No. These are replacement components for medical training simulators only and must not be used for any human dental application.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.