



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Stomach

**Item No.** LF03606U  
**Weight** 0.02kg  
**Make:** Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:** Airway Management Adult

## Product Description

For use with **Life/form®** Airway Management Trainers and **Life/form®** Adult **CRiSis™** Manikins.

### Why Choose the Life/form® Stomach for Your Simulation Lab in India?

The Life/form® Stomach (Item No. LF03606U / W44054) by Nasco Healthcare, USA is a 0.02 kg genuine OEM replacement stomach bladder for Life/form® Airway Management Trainers and Life/form® Adult CRiSis™ Manikins — compatible with W44104 ("Airway Larry" Airway Management Trainer with Stand / LF03699U), W44124 (Complete Adult CRiSis™ Resuscitation Training System Full-Body), and W44135 (Deluxe CRiSis™ Torso Manikin). The stomach bladder is the consumable component that inflates visibly during oesophageal intubation or excessive BVM ventilation pressure — the direct consequence-feedback mechanism that teaches trainees to verify tube placement before

proceeding and to calibrate BVM pressure correctly.

The stomach bladder's training function is specific and irreplaceable: it makes the consequence of oesophageal intubation immediately visible — the stomach swells, the lungs do not inflate — providing the real-time feedback that verbal instruction alone cannot achieve at the moment of error. Without an intact stomach bladder in place, the oesophageal intubation consequence feedback is absent and the trainer cannot deliver its full clinical training value. The bladder degrades from repeated inflation cycles and must be replaced when visible distension becomes inconsistent or the bladder fails to hold inflation. The 0.02 kg weight reflects the minimal mass of this thin inflatable consumable — easy to stock in multiples without logistics overhead. Only genuine OEM LF03606U material is guaranteed to fit the compatible trainer airway modules with correct connection geometry and correct inflation characteristics.

**Why Buy from SEM Trainers?** SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

### Frequently Asked Questions

**Q: Which Life/form® trainers is this stomach compatible with?** A: LF03606U is compatible with W44104 ("Airway Larry" Airway Management Trainer with Stand / LF03699U), W44124 (Complete Adult CRiSis™ Full-Body), and W44135 (Deluxe CRiSis™ Torso Manikin) by Nasco Healthcare. Confirm compatibility with your specific trainer configuration before ordering.

**Q: Why does the stomach bladder need to be replaced?** A: Repeated inflation cycles from oesophageal intubation training degrade the bladder material, causing inconsistent inflation or failure to hold distension. A degraded stomach bladder removes the oesophageal intubation consequence feedback that is central to this trainer's clinical training value.

**Q: How do I order and get support from SEM Trainers?** A: Contact us at [sem@semtrainers.com](mailto:sem@semtrainers.com) or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.