



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE


SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Lungs

Item No. LF03605U

Weight 0.08kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Airway Management Adult

Product Description

for W44104, W44124, W44134, W44135

The Life/form® Airway Larry Lungs (Item No. LF03605U / 1005607) by Nasco Healthcare, USA is a 0.08 kg genuine OEM replacement lung set for Life/form® Airway Larry Airway Management Trainers, compatible with W44104, W44124, W44134, and W44135. The bifurcated lungs in Life/form® Airway Larry trainers serve two critical training functions simultaneously: they inflate visibly with successful ventilation, confirming correct tracheal tube placement, and they enable bilateral auscultation that distinguishes correct bilateral lung expansion from right mainstem or oesophageal intubation. When the lungs degrade — losing their correct compliance, developing micro-perforations, or failing to inflate symmetrically — both of these feedback mechanisms are compromised, and trainees receive

inaccurate confirmation of tube placement.

For simulation labs running regular intubation training sessions across multiple trainees, lung replacement is the most common maintenance task for Airway Larry trainers. The 0.08 kg weight reflects the compact, low-mass nature of the lung consumable — easy to stock in quantity without logistics overhead. Only genuine OEM LF03605U lungs are guaranteed to fit the compatible Airway Larry trainer airway modules with correct dimensions, correct compliance characteristics, and correct inflation response. Using non-OEM substitutes risks incorrect inflation resistance, asymmetric lung compliance, and poor bilateral auscultation clarity — all of which degrade the quality of tube placement feedback that makes the bifurcated lung design clinically valuable in the first place.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which Airway Larry trainers are these lungs compatible with? A: LF03605U is compatible with Life/form® Airway Larry Airway Management Trainers W44104 (LF03699U), W44124, W44134, and W44135 by Nasco Healthcare. For other Life/form® trainers, see the Replacement Lungs/Stomach (W99837).

Q: How do degraded lungs affect training quality? A: Degraded lungs cause incorrect ventilation confirmation feedback — they may fail to inflate visibly on correct intubation or inflate asymmetrically, giving trainees false confirmation or false negative results. Bilateral auscultation clarity also degrades. Replace when inflation consistency is reduced or visible asymmetry occurs.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.