



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Cricothyrotomy Simulator

Item No. LF00994
Weight 6.7kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Advanced Trauma Life Support (ATLS)

Product Description

The newly designed **Life/form®** Cricothyrotomy Simulator has been developed for learning and practicing the techniques necessary to perform needle or surgical cricothyrotomy procedures. Paramedics, EMTs, combat medics, flight nurses, anesthesiologists, and other emergency medical personnel will have the opportunity to strengthen their ability and confidence to perform or assist in implementing surgical airways.

- Anatomically accurate landmarks aid in site training and allow for fast action.
- The hyperextended neck allows the user to determine the proper incision site. The trachea in the simulator is replaceable, as the airway passes completely through from top to bottom. This allows checking the stylet and obturator placement once the incision has been made.
- Complete with a chin and full-size neck, ties can be used to hold the obturator in a secure position. Inflation of the simulated lung verifies correct placement.

Includes the simulator with base, six replaceable neck skins, six adult trachea inserts (four rigid, two soft), six child trachea inserts (four rigid, two soft), two simulated lungs, instruction manual, and hard carry case.

	Deluxe Cricothyrotomy Simulator	TruCric	Cricotracheotomy Trainer	Cricothyrotomy Simulator
Item No.	1017948	1021893	1005180	1020255
Surgical cricothyrotomy	Yes	Yes	Yes	Yes
Needle cricothyrotomy	Yes	Yes	Yes	Yes
Insertion of endotracheal tube	Yes	Yes	Yes	Yes
Full use of devices such as the PCK and UPK	-	Yes	Yes	Yes
Seldinger technique	-	Yes	Yes	Yes
Anatomical Landmarks	Cricothyroid membrane	Cricothyroid membrane, laryngeal cartilages and tracheal rings	Cricothyroid membrane, chin and full-sized neck	Cricothyroid membrane, chin and full-sized neck
Review airway anatomy	-	-	Yes	Yes
Realistic soft skin	Yes	Yes	-	-
Confirm proper placement of endotracheal tube	Inflation of lung indicates proper placement	Inflation of lung indicates proper placement	no	Inflation of lung indicates proper placement
Simulation features	• Allows for an adult as well as pediatric cricothyrotomy • Includes adult and 3-year-old tracheas	• Neck skin can be rotated on the model allowing up to 20 incisions • Larynx insert allows 1 surgical cricothyroidotomy and 3 tracheostomy procedures	• Up to 36 incisions on each unit • Up to 18 mini tracheostomies	• Assessment of styled and obturator placement (replaceable trachea)