



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Amputated Bleeding Arm for SMART STAT

Item No. 1018592 [W44181]
Weight 0.342kg
MPN: 101-468U / PP00468U / SB48724U
Previous Item No. 3004290

[Read More](#)

SKU:

Categories:ALS Adult



Product Description

SMART STAT has the capability to put in place this amputation mouldage so practitioners have the opportunity to treat this trauma with a tourniquet (TQ). When the TQ is applied correctly, the bleeding will cease. Blow flow is controlled by the instructor by selecting the function on the PPC (Included with SMART STAT)

Why Choose the Amputated Bleeding Arm for SMART STAT for Your Simulation Lab in India?

The Amputated Bleeding Arm for SMART STAT (Item No. 1018592, MPN: 101-468U) by Nasco Healthcare, USA is a purpose-built amputation mouldage accessory designed exclusively for the Simulaids SMART STAT ALS simulator. Weighing just 0.342 kg, it attaches to the SMART STAT to create a realistic traumatic arm amputation scenario — complete with active bleeding that trainees

must control using a tourniquet. This is one of the most clinically relevant trauma training scenarios in pre-hospital and emergency medicine education, and this accessory makes it possible to train it repeatedly, consistently, and with meaningful physical feedback on the SMART STAT platform your lab already owns.

The mechanism is straightforward and instructor-controlled: blood flow through the amputated arm is activated and managed by the instructor via the SMART STAT's PPC (Portable Patient Controller), which is included with the SMART STAT simulator. When a trainee applies a tourniquet correctly — with proper placement and adequate tension — the bleeding stops. This creates an objective, measurable training outcome that reinforces correct tourniquet technique and confirms competency in haemorrhage control, a core skill requirement in TCCC (Tactical Combat Casualty Care), ATLS, and pre-hospital trauma training programmes across India. The accessory is durable, easy to attach and remove, and integrates directly with the SMART STAT's existing bleeding system without requiring additional hardware or complex setup.

This accessory is ideal for ATLS and TCCC training programmes, pre-hospital and paramedic training institutes, armed forces and paramilitary medical training units, emergency medicine departments, and disaster response teams across India that use the SMART STAT simulator and want to add high-realism tourniquet and haemorrhage control training to their curriculum.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare and Simulaids products in India with over 30 years of experience supplying medical simulators and accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Is this accessory compatible with all SMART STAT configurations? A: Yes. This amputated bleeding arm is designed for use with the Simulaids SMART STAT ALS simulator. Bleeding is controlled via the SMART STAT's PPC (Portable Patient Controller), which is included with the SMART STAT unit.

Q: Does tourniquet application actually stop the bleeding during training? A: Yes. When the trainee applies a tourniquet correctly, the bleeding from the amputated arm ceases — providing immediate, objective feedback on technique.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare accessories with pan-India delivery.