



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE


SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Defibrillation Chest Skin

Item No. LF03611U

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Adult

Product Description

Perfect for practicing defibrillation on any Child CRiSis™ or Child Airway Management Torso. The chest skin lets you defibrillate using standard manual, automatic, or semi-automatic external monitor defibrillators by attaching directly to your manikin. The skin also makes it possible to monitor the manikin, just like a real patient.

Why Choose the Life/form® Defibrillation Chest Skin for Your Simulation Lab in India?

The Life/form® Defibrillation Chest Skin (Item No. LF03611U) by Nasco Healthcare, USA is the genuine OEM defibrillation upgrade skin designed for the Child CRiSis™ manikin and Child Airway Management Torso. It adds full defibrillation and patient monitoring capability to these child manikin platforms — allowing trainees to practice paediatric defibrillation using real manual, automatic, and

semi-automatic external defibrillators attached directly to the manikin's chest. This is a critical training capability for any simulation lab that runs PALS (Paediatric Advanced Life Support) or paediatric emergency medicine courses, where hands-on practice with real defibrillation equipment on an appropriately sized manikin makes the difference between a credible scenario and a compromised one.

The chest skin works by converting any compatible Child CRiSis™ or Child Airway Management Torso into a defibrillation-capable trainer. Real defibrillator pads or paddles attach directly to the conductive chest skin, allowing the delivery of actual defibrillation energy — not simulated. This means trainees develop pad placement habits, charging and delivery procedures, and rhythm recognition skills with the same equipment they will use in clinical practice. The skin also enables patient monitoring — trainees can connect the manikin to a cardiac monitor exactly as they would a real paediatric patient, integrating monitoring and defibrillation decision-making into a complete PALS simulation scenario.

This skin is ideal for PALS training programmes, paediatric emergency medicine departments, paediatric nursing and medical skills labs, and NMC-mandated simulation labs in Indian medical colleges that train paediatric resuscitation competencies.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which manikins is this defibrillation chest skin compatible with? A: This chest skin (LF03611U) is designed for the Child CRiSis™ manikin and the Child Airway Management Torso by Nasco Healthcare. It is not compatible with adult CRiSis™ configurations.

Q: Can real defibrillator equipment be used with this skin? A: Yes. Standard manual, automatic, and semi-automatic external defibrillators attach directly to the skin, and the manikin can be monitored like a real patient.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare accessories with pan-India delivery.