

**SEM****SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724**Mobile 1:** +91 98791 03905**Email:** sem@semtrainers.com**Website:** www.semtrainers.com

Defibrillation Pad & Patient Adapter Physio Control Adapters for CRiSis Manikin

Item No. LF03656U**Make:** Nasco Healthcare, USA[Read More](#)**SKU:****Categories:** ALS Adult

Product Description

The reusable defibrillation pad adapters are designed to be used in place of expensive defib gel pads. They come complete with patient simulator adapters that will enable you to use the patient simulator you already have with the CRiSis™ Manikin, CRiSis™ Manikin Update Package, or Defibrillation Chest Skin. The training cables to use with the CRiSis™ manikin are available with unique cable ends for use with the following defibrillators: Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick combo. If your defibrillator uses snap-on patches we also have the adapters you need to eliminate the use of pads. Training cables are not recommended for use with non-rechargeable battery AED. This includes FirstSave, HeartStream, and Medtronic Physio control with non-rechargeable lithium batteries. This package will save you both time and money!

Why Choose the Defibrillation Pad & Patient Adapter (Physio Control Adapters) for CRiSis™

Manikin for Your Simulation Lab in India?

The Defibrillation Pad & Patient Adapter for CRiSis™ Manikin (Item No. LF03656U) by Nasco Healthcare, USA is a reusable adapter set that eliminates the recurring cost of disposable defibrillation gel pads during CRiSis™ manikin training sessions. For simulation labs running regular ACLS and ALS training programmes across India, disposable gel pads represent a significant and avoidable ongoing expense. This package replaces them with reusable pad adapters that connect directly to the CRiSis™ manikin's defibrillation chest skin — delivering the same training functionality at a fraction of the per-session consumable cost.

The package includes both pad adapters and patient simulator adapters, enabling the CRiSis™ manikin to interface with the patient simulator already in use in your lab. The training cables supplied cover the four defibrillator brands most commonly found in Indian hospital simulation labs and training centres: Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick Combo. For defibrillators that use snap-on patches rather than gel pads, snap-on adapters are also included to eliminate the need for pads entirely. This package is compatible with the CRiSis™ Manikin, CRiSis™ Manikin Update Package, and the Defibrillation Chest Skin. One important note: these training cables are not recommended for use with non-rechargeable battery AEDs, including FirstSave, HeartStream, and Medtronic Physio Control units with non-rechargeable lithium batteries.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which defibrillator brands are the training cables compatible with? A: The package includes cable ends compatible with Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick Combo defibrillators. Snap-on patch adapters are also included for defibrillators that use snap-on patches.

Q: Are these adapters compatible with all AED models? A: No. These training cables are not recommended for use with non-rechargeable battery AEDs, including FirstSave, HeartStream, and Medtronic Physio Control units with non-rechargeable lithium batteries.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare accessories with pan-India delivery.