



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Life/form® CRiSis™ Manikin Replacement Lungs

Item No. LF03285U

Weight 0.114kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Adult

Product Description

For use with “Airway Larry”, Adult CRiSis™ and CPARLENE® Manikins

Life/form® CRiSis™ Manikin Replacement Lungs (Item No. LF03285U | Weight: 0.114 kg | Nasco Healthcare, USA)

Genuine OEM replacement lung bags for the Life/form® CRiSis™ family of airway training manikins. Lung bags are high-use consumables — they undergo repeated compression during ventilation practice and will eventually develop micro-tears, loss of elasticity, or seal failure that degrades the realism of breathing simulation. Replacing them promptly restores correct lung compliance, accurate chest rise, and consistent ventilation resistance to original factory performance.

Compatible with Airway Larry®, Adult CRiSis™, and CPARLENE® manikins. [semtrainers](https://semtrainers.com)

Key Features:

- Genuine Life/form® OEM replacement part — ensures correct fit, compliance, and seal integrity
- Restores realistic chest rise and ventilation resistance for high-fidelity airway management training
- Lightweight at 0.114 kg — easy to store and swap during training sessions
- Compatible across the CRiSis™ manikin family: Airway Larry®, Adult CRiSis™, and CPARLENE®

When to replace:

- Visible cracking, deformation, or loss of shape in the lung bag material
- Inconsistent or reduced chest rise during bag-mask or mouth-to-mouth ventilation
- Air leakage detected during ventilation assessment
- As part of scheduled preventive maintenance — recommended after high-frequency use periods

Recommended practice: Keep at least one spare set on hand. Lung bag failure mid-session halts training. With a replacement readily available, downtime is measured in minutes, not days waiting for shipment.

For bulk consumable orders or multi-unit institutional procurement, contact SEM Trainers.