



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Lung Bag Connector & O-Ring for 1005650

Item No. 1019770
MPN: LF03321U

[Read More](#)

SKU:

Categories:ALS Adult

Product Description

Lung Bag Connector & O-Ring for Manikin 1005650 (Item No. 1019770 | MPN: LF03321U)

A genuine OEM replacement connector and O-ring set designed specifically for the lung bag assembly used in manikin model 1005650. When the lung bag connector develops wear or the O-ring loses its seal, ventilation resistance changes and breathing simulation fidelity degrades — this kit restores correct airway function and airtight lung bag operation to original factory specification.

Fits: Manikin model 1005650 **Includes:** Lung bag connector + O-ring **Purpose:** Restores proper lung bag seating and airtight ventilation simulation **Recommended use:** Replace when ventilation resistance feels inconsistent, when visible O-ring cracking or deformation is present, or as part of

scheduled preventive maintenance

Keeping spare connectors and O-rings on hand avoids training downtime. Order before failure, not after.

For compatibility queries or bulk spare parts requirements, contact SEM Trainers.