



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Trunk shell for CRiSis™ ALS simulator

Item No. 1019762

MPN: LF03190U

[Read More](#)

SKU:

Categories:ALS Adult

Product Description

The trunk shell is used as a replacement part for the ALS (advanced life support) simulator CRiSis™ W44135.

Why Choose the Life/form® Trunk Shell for CRiSis™ ALS Simulator for Your Simulation Lab in India?

The Life/form® Trunk Shell for CRiSis™ ALS Simulator (Item No. 1019762, MPN: LF03190U) by Nasco Healthcare, USA is the genuine OEM replacement trunk shell for the CRiSis™ ALS simulator family. The trunk shell is the structural outer body of the CRiSis™ torso — the component that houses and protects the internal CPR mechanism, chest decompression modules, defibrillation skin mount points, and IV and BP arm attachment interfaces. When the trunk shell sustains damage from high-

frequency training use, impact, or material fatigue, replacing it with the genuine OEM part is the correct approach to restore the manikin to full structural and functional specification.

Attempting to continue training on a damaged or cracked trunk shell is not simply a cosmetic issue. Structural integrity of the trunk shell affects CPR compression feedback — a deformed shell can alter the resistance and recoil feel that trainees calibrate their technique against. It also affects the seating and alignment of the defibrillation chest skin, IV arm, and BP arm modules, any of which can become misaligned or functionally compromised if the shell is distorted. For simulation labs running high-throughput BLS, ACLS, or ALS certification programmes — where the CRiSis™ manikin may undergo dozens of compression cycles and procedural interventions per training day — the trunk shell is a structural consumable with a finite service life. Maintaining a spare LF03190U in stock prevents unplanned training downtime when shell replacement becomes necessary. Only the genuine OEM trunk shell guarantees correct dimensional fit, correct mounting geometry for all CRiSis™ modules, and material compliance with the original manikin specification.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM spare parts to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which CRiSis™ simulators is this trunk shell compatible with? A: This trunk shell (LF03190U, Item No. 1019762) is designed as a replacement part for the Life/form® CRiSis™ ALS simulator family by Nasco Healthcare. Contact us to confirm compatibility with your specific CRiSis™ configuration before ordering.

Q: Does replacing the trunk shell require any special tools or technician support? A: SEM Trainers provides on-site technical support for CRiSis™ manikin servicing across India. Contact us to arrange support if needed: sem@semtrainers.com or +91-88495 63724.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare spare parts with pan-India delivery.