



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Deluxe CRiSis™ Torso Manikin

Item No. LF03958U
Weight 10.493kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Adult

Product Description

Life/form® Deluxe CRiSis™ Torso Manikin features the same quality, durability, and reliability as the Deluxe CRiSis™ Manikin (W44128) only in a smaller, easier-to-transport version. Unlike full-body manikins, this ACLS torso is much smaller and easier to move from facility to facility for training. This system now features the “Airway Larry” Airway Management Training Head with carotid pulse, Defibrillation Chest Skin, and Interactive ECG Simulator (W44153), all attached to a CPARLENE® Torso. Includes a 8-oz. pump spray bottle of lubricant (W44105) and comes in a durable, foam-lined carton (hard carry case is available, (W44136).

Why Choose the Life/form® Deluxe CRiSis™ Torso Manikin for Your Simulation Lab in India?

The Life/form® Deluxe CRiSis™ Torso Manikin (Item No. LF03958U) by Nasco Healthcare, USA is the compact, ECG-enabled torso configuration of the CRiSis™ platform — and the key differentiator

between this and the standard CRiSis™ Torso is the inclusion of the Interactive ECG Simulator as standard. Weighing 10.493 kg, it integrates the "Airway Larry" Airway Management Training Head with carotid pulse, Defibrillation Chest Skin, and Interactive ECG Simulator — all mounted on a CPARLENE® Torso — in a single portable unit that is explicitly designed to be easier to transport between facilities than a full-body manikin. It ships with a lubricant spray bottle and a foam-lined protective carton; a hard carry case (LF03467U) is available separately.

The "Airway Larry" head provides anatomically accurate airway management training — oral, nasal, and digital intubation with ET, EOA, PTL, Combitube®, King Systems, and supraglottic airway devices, suction technique, palpable cricoid cartilage for Sellick manoeuvre, and palpable manual carotid pulse. The Defibrillation Chest Skin is compatible with all standard defibrillator brands and types, with ECG monitoring capability for realistic patient monitoring scenarios. The Interactive ECG Simulator adds programmed rhythm scenarios for rhythm recognition and defibrillation decision training integrated into ALS scenarios — the capability that separates this torso from the standard CRiSis™ Torso (LF03957U). An IV arm (available separately) can be added to expand the training scope to IV access when needed.

This torso is the right choice for ALS training programmes, emergency medicine departments, outreach instructor programmes, and NMC-mandated simulation labs across India that need a highly portable, ECG-capable airway and resuscitation trainer that can move between campuses or training sites without the logistical burden of a full-body system.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: How does this differ from the CRiSis™ Torso (LF03957U)? A: The Deluxe CRiSis™ Torso includes the Interactive ECG Simulator as standard — the standard CRiSis™ Torso does not. Both include the "Airway Larry" head and Defibrillation Chest Skin. If ECG rhythm training is part of your ALS curriculum, this is the correct choice.

Q: Can an IV arm be added to this torso? A: Yes. The IV Arm for 1005650 is available as a separate accessory and attaches to this torso to add IV access training capability.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.

Manikin can be used with Supraglottic airway devices for airway training.

