



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



IV Arm Replacement Skin/Veins

Item No. LF03641U

Weight 0.06kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Child

Product Description

For use with Life/form® Infant CRiSis™ Manikins. Caution: Use only Life/form® Replacements. DO not use conventional latex tubing. This will cause malfunctioning of the simulator.

Why Choose the Life/form® IV Arm Replacement Skin/Veins for Infant CRiSis™ Manikins for Your Simulation Lab in India?

The Life/form® IV Arm Replacement Skin/Veins (Item No. LF03641U) by Nasco Healthcare, USA is the genuine OEM combined replacement kit — skin and veins together — for the IV arm module on Life/form® Infant CRiSis™ Manikins, weighing 0.06 kg. The Infant CRiSis™ IV arm features smaller lumen tubing that simulates the size of an infant's veins, providing access to the cephalic and basilic veins and the dorsal venous arch on the hand with realistic pressurised flashback. This combined kit

restores both consumable components of the IV arm simultaneously — the outer skin layer that provides surface resistance and cannulation entry, and the internal vein tubing that provides vein palpability, correct flashback behaviour, and intraluminal resistance — to full factory specification after high-frequency training use.

The skin and veins on the Infant CRiSis™ IV arm are proportionally smaller than those on child and adult IV training arms, reflecting the significantly finer vessel calibre that paediatric venipuncture requires. This smaller lumen tubing is not interchangeable with replacement kits from other manikin families — only genuine OEM LF03641U material is guaranteed to match the correct infant lumen dimensions, correct flashback pressure, and correct material compliance for the Infant CRiSis™ IV arm. Using a non-OEM substitute risks incorrect vein geometry, inconsistent flashback, and simulator malfunction. **Nasco Healthcare explicitly cautions that only genuine Life/form® OEM replacement components must be used with Infant CRiSis™ Manikins — conventional latex tubing will cause malfunction.** For simulation labs running regular neonatal and paediatric IV access training sessions on the Infant CRiSis™, maintaining a stock of LF03641U prevents unplanned training interruption and ensures consistent training quality for trainees developing the fine-motor IV access skills that infant venipuncture demands.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which Infant CRiSis™ simulators is this replacement skin/vein kit compatible with? A: LF03641U is designed for the IV arm module on Life/form® Infant CRiSis™ Manikins by Nasco Healthcare. Contact us to confirm compatibility with your specific configuration before ordering.

Q: Does this kit include both the skin and the veins? A: Yes. LF03641U is a combined kit containing both the replacement outer skin and the replacement internal vein tubing for the Infant CRiSis™ IV arm.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.