

**SEM****SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724**Mobile 1:** +91 98791 03905**Email:** sem@semtrainers.com**Website:** www.semtrainers.com

Femoral Injection Pad

Item No. LF03626U**Weight** 0.1kg**Make:** Nasco Healthcare, USA[Read More](#)**SKU:****Categories:**ALS Child

Product Description

An injection simulator for W44091 and W44126. The femoral injection pad is a great alternative to learning on real patients, allowing realistic and repeated injections.

Why Choose the Life/form® Femoral Injection Pad for Your Simulation Lab in India?

The Life/form® Femoral Injection Pad (Item No. LF03626U) by Nasco Healthcare, USA is a genuine OEM replacement injection pad for the femoral access module on compatible Life/form® Child CRiSis™ manikins, weighing 0.1 kg. It is compatible with the Deluxe Child CRiSis™ Manikin with ECG Simulator (W44091) and W44126. The femoral injection pad maintains the realistic injection simulation capability of the Child CRiSis™ femoral access leg — the module that provides trainees with palpable symphysis pubis and anterior superior iliac spine landmarks, palpable arterial pulse, and pressurised

venous flashback for femoral vascular access training in paediatric PALS scenarios.

Femoral venous access is a critical paediatric emergency skill when peripheral and IO access have failed or are contraindicated, and correct technique depends on accurate landmark identification and correct angle of insertion relative to the femoral artery — a high-stakes skill where training on a worn or degraded pad translates directly into reduced landmark accuracy and incorrect tactile feedback. The femoral injection pad provides the tissue resistance and surface presentation that trainees calibrate their technique against during each session. As with all injection site consumables, repeated needle insertions progressively degrade the pad material, reducing its ability to accurately represent the tissue resistance and surface anatomy of the femoral triangle. Replacing the pad regularly with genuine OEM LF03626U material restores the correct surface feel, correct landmark presentation, and correct injection resistance that paediatric femoral access training requires. Only genuine OEM replacement pads are guaranteed to fit the femoral access module of the compatible Child CRiSis™ manikins with correct dimensions and material specification.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which Child CRiSis™ manikins is this femoral injection pad compatible with? A: LF03626U is compatible with the Life/form® Deluxe Child CRiSis™ Manikin with ECG Simulator (W44091) and W44126 by Nasco Healthcare. Contact us to confirm compatibility with your specific configuration before ordering.

Q: How does the femoral injection pad degrade and when should it be replaced? A: Repeated needle insertions progressively degrade the pad material, reducing tissue resistance accuracy and surface landmark definition. Replace when the pad shows visible puncture accumulation, loss of elasticity, or no longer provides consistent resistance feedback during injection training.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.