



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Life/form® Injectable Training Arm Replacement Skin and Vein Kit

Item No. LF03629U

Weight 0.42kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Child

Product Description

For use with the **Life/form®** Child CPR/Airway Management Torso (LF03633).

Why Choose the Life/form® Injectable Training Arm Replacement Skin and Vein Kit for Your Simulation Lab in India?

The Life/form® Injectable Training Arm Replacement Skin and Vein Kit (Item No. LF03629U) by Nasco Healthcare, USA is the genuine OEM combined replacement kit — skin and veins together — weighing 0.42 kg, designed specifically for use with the Life/form® Child CPR/Airway Management Torso (LF03633). It restores the injectable arm module on that manikin to full factory specification, replacing both the outer skin and the internal vein tubing simultaneously after high-frequency IV access training use. This kit is not interchangeable with the adult injectable training arm replacement kit (LF03215U) —

the child torso's arm has different dimensions, different vein channel geometry, and smaller lumen tubing that reflects the proportionally finer vessel calibre of a paediatric patient.

The injectable arm on the Life/form® Child CPR/Airway Management Torso is a dedicated paediatric IV access training module. The outer skin provides the surface tissue resistance and surface landmark feel that trainees must work against during cannulation; the internal vein tubing provides the correct lumen size, correct flashback behaviour, and correct intraluminal resistance that a child-scale IV training arm requires. Both components degrade from repeated needle insertions — the skin accumulating puncture damage that destroys surface realism, the veins developing micro-perforations that impair flashback consistency. Replacing both together with genuine OEM LF03629U material restores the correct paediatric IV access training experience. Only genuine OEM material is guaranteed to fit the LF03633 arm module with correct dimensions, correct vein geometry, and correct material compliance. For simulation labs running regular paediatric nursing, PALS, or emergency medicine IV access training sessions on the Child CPR/Airway Management Torso, maintaining a stock of LF03629U prevents training interruption and protects training quality.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which manikin is this kit compatible with? A: LF03629U is designed for the injectable arm module on the Life/form® Child CPR/Airway Management Torso (LF03633) by Nasco Healthcare. It is not compatible with the adult injectable training arm models.

Q: Is this the same kit as the adult Injectable Training Arm Replacement Skin and Vein Kit (LF03215U)? A: No. LF03629U is specifically sized for the child torso (LF03633) with smaller lumen tubing reflecting paediatric vein calibre. LF03215U is for adult injectable training arm models. They are not interchangeable.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.