



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



2 Intraosseous Leg Skins

Item No. LF03619U

Weight 0.22kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Child

Product Description

Package of 2.

Why Choose the Life/form® 2 Intraosseous Leg Skins for Your Simulation Lab in India?

The Life/form® 2 Intraosseous Leg Skins (Item No. LF03619U) by Nasco Healthcare, USA is a 2-pack of genuine OEM replacement outer leg skins for the intraosseous access module on compatible Life/form® paediatric ALS simulators, weighing 0.22 kg per pack. The IO leg skin is the external tissue layer through which the IO needle or drill first passes during intraosseous access training — it is what the trainee's hands contact first, and it provides the initial tissue resistance and surface landmark feel that correct IO technique depends on. It is a distinct consumable from the IO Leg Bones (LF03618U), which are the inner cortical bone inserts that provide the deeper resistance of penetrating simulated

periosteum and cortical bone.

The skin and bone components of an IO training module degrade in different ways and at different rates. The leg skin accumulates puncture damage from repeated IO needle insertions at the skin surface — the tissue layer loses its integrity, develops visible puncture tracks, and eventually fails to provide the correct initial resistance and surface texture that trainees calibrate their technique against. A worn skin also misrepresents the anatomical surface landmarks — the tibial tuberosity, flat medial surface, and correct insertion point geometry — that trainees must learn to identify reliably for safe IO placement in real paediatric patients. Replacing the leg skins with genuine OEM LF03619U skins restores correct surface resistance, correct anatomical landmark presentation, and correct tissue feel. For simulation labs running regular PALS, paediatric emergency medicine, and NLS training programmes, both the skin (LF03619U) and bone (LF03618U) consumables should be stocked and replaced on independent schedules — the skin typically degrades faster than the bone at the surface entry point.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which Life/form® simulators are these IO leg skins compatible with? A: LF03619U is designed for compatible Life/form® paediatric ALS simulators with intraosseous tibial access modules by Nasco Healthcare. Contact us to confirm compatibility with your specific simulator before ordering.

Q: What is the difference between this skin pack (LF03619U) and the IO Leg Bones (LF03618U)? A: The leg skin (LF03619U) is the outer tissue layer providing surface resistance and anatomical landmark presentation. The leg bone (LF03618U) is the inner cortical insert providing deeper bone penetration resistance. Both are separate consumables that degrade at different rates and must be stocked and replaced independently.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.