



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Adapter with Control Training Cables

Item No. 1017990  
Weight 0.05kg  
MPN: LF03962U

[Read More](#)

**SKU:**

**Categories:**BLS and CPR Accessories

## Product Description

The reusable defibrillation pad adapters are designed to be used in place of expensive defib gel pads. They come complete with patient simulator adapters that will enable you to use the patient simulator you already have with the CRiSis™ Manikin, CRiSis™ Manikin Update Package, or Defibrillation Chest Skin. The training cables to use with the CRiSis™ manikin are available with unique cable ends for use with the following defibrillators: Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick combo. If your defibrillator uses snap-on patches we also have the adapters you need to eliminate the use of pads. Training cables are not recommended for use with non-rechargeable battery AED. This includes FirstSave, HeartStream, and Medtronic Physio control with non-rechargeable lithium batteries. This

package will save you both time and money!

## Why Choose the Adapter with Control Training Cables for Your Simulation Lab in India?

The Adapter with Control Training Cables (Item No. 1017990, MPN: LF03962U) by Nasco Healthcare, USA is a 0.05 kg reusable defibrillation training accessory for the Life/form® CRiSis™ Manikin, CRiSis™ Manikin Update Package, or Defibrillation Chest Skin — designed to replace expensive single-use defib gel pads with reusable pad adapters that connect your existing defibrillator to the CRiSis™ platform without consumable pads. Patient simulator adapters are included, enabling the defibrillator you already use in training to work directly with the CRiSis™ system. For simulation labs, ACLS training centres, and hospital in-service programmes in India running regular defibrillation training sessions on the CRiSis™ platform, this adapter set eliminates the recurring cost and logistics of sourcing and stocking single-use gel pads — delivering cost savings and training continuity across sessions.

Training cables are supplied with unique cable ends matched to specific defibrillator models — Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick Combo — and adapters for defibrillators using snap-on patches are also available, eliminating pad use entirely on those platforms.

**Important compatibility restriction:** Training cables are not recommended for use with non-rechargeable battery AEDs. This includes FirstSave, HeartStream, and Medtronic Physio Control units with non-rechargeable lithium batteries — confirm your defibrillator model and battery type before ordering. Contact SEM Trainers if you are unsure whether your defibrillator is compatible.

**Why Buy from SEM Trainers?** SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

## Frequently Asked Questions

**Q: Which manikin systems is this adapter compatible with?** A: This adapter set is designed for the Life/form® CRiSis™ Manikin, CRiSis™ Manikin Update Package, and Defibrillation Chest Skin by Nasco Healthcare.

**Q: Which defibrillators are supported?** A: Training cables are available for Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick Combo. Adapters for snap-on patch defibrillators are also available. Not recommended for non-rechargeable battery AEDs including FirstSave, HeartStream, or Medtronic Physio Control with non-rechargeable lithium batteries.

**Q: How do I order and get support from SEM Trainers?** A: Contact us at [sem@semtrainers.com](mailto:sem@semtrainers.com) or call +91-88495 63724 to confirm compatibility with your defibrillator before ordering.