



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



CPR Simon Full Body Simulator w/ OMNI® Code Blue Pack

Item No. 1009220 [W45116]
Weight 20.41kg
MPN: S311.M

[Read More](#)

SKU:

Categories:BLS Adult



Product Description

Quality CPR simulator designed to operate for years! Strong, sturdy construction with full-size, realistic landmarks and one piece disposable airways. Prewired for CPRLink™ for a first class BLS program. Features - Totally disposable airways to practice hygienic mouth-to-mouth ventilation - Fully articulating

head, neck, and jaw permitting head/tilt, chin lift, jaw thrust, and neck extension into the “sniffing” position - Airway blocked when head is forward - Easily accessible chest cavity with ribcage, lungs, and heart for realism and practice of the Heimlich maneuver - Realistic chest rise during ventilation - Carotid, femoral arterial pulse sites - Femoral venous site · Intramuscular (IM) Injection Sites on the deltoids, quadriceps and left gluteus medius - Eyes open and close; one pupil is dilated - Omni® Code Blue® pack monitors and logs the cadence and depth of cardiac compression and airway ventilation - Jointed elbows, wrists, knees and ankles - 10 disposable airways - Detachable at waist for easy storage - Instruction Manual - Carrying bag

Why Choose the CPR Simon Full Body Simulator with OMNI® Code Blue Pack for Your Simulation Lab in India?

The CPR Simon Full Body Simulator with OMNI® Code Blue Pack (Item No. 1009220 / W45116, MPN: S311.M) by Gaumard Scientific, USA is a 20.41 kg, strong, sturdy full-body adult simulator built to operate for years, with full-size realistic anatomical landmarks and one-piece disposable airways. Pre-wired for CPRLink™, this simulator is designed as the foundation of a first-class BLS programme — for ACLS training centres, hospital in-service programmes, and nursing institutes in India that need a comprehensive full-body platform spanning CPR, airway management, vascular access, and injection technique on a single durable manikin.

The fully articulating head, neck, and jaw permit head tilt/chin lift, jaw thrust, and neck extension into the “sniffing” position for advanced airway training — and the airway blocks correctly when the head is positioned forward, training trainees to recognise incorrect positioning. An easily accessible chest cavity with ribcage, lungs, and heart adds realism and supports practice of the Heimlich manoeuvre. Realistic chest rise confirms effective ventilation. Carotid and femoral arterial pulse sites support pulse assessment training, a femoral venous site enables vascular access practice, and intramuscular injection sites on the deltoids, quadriceps, and left gluteus medius round out the clinical procedure training capability. Responsive eyes open and close, with one pupil dilated for neurological assessment.

The OMNI® Code Blue® pack monitors and logs the cadence and depth of cardiac compression and airway ventilation, giving instructors objective CPR quality data during training. Jointed elbows, wrists, knees, and ankles allow realistic patient positioning and movement during scenario-based training. The manikin detaches at the waist for easier storage and transport. Includes 10 disposable airways, an instruction manual, and a carrying bag.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Gaumard Scientific products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Frequently Asked Questions

Q: What does the OMNI® Code Blue Pack monitor? A: It monitors and logs the cadence and depth of cardiac compression and airway ventilation, giving instructors objective CPR quality data during training sessions.

Q: What vascular access and injection sites are included? A: A femoral venous site for vascular access practice, plus intramuscular injection sites on the deltoids, quadriceps, and left gluteus medius.

Q: Can the manikin be disassembled for storage? A: Yes. It detaches at the waist for easier storage and transport, and is supplied with a carrying bag.