



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Baby Buddy™ Single CPR Manikin

Item No. LF03720

Weight 0.4kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:BLS Newborn

Product Description

• Full-Body Infant Manikin with Proper Landmarks to Realistically Practice CPR Techniques • Sanitary One-Piece Lung/Mouth Protection System • Economical - Each Student Can Have Their Own Manikin • Low-Cost Disposable Lung/Mouth Bags - No Messy Cleanup • Lung/Mouth Bag Installs Easily (Tool Included) • Removable Chest for Easy Lung/Mouth Bag Installation • Visible Chest Rise When Ventilated • Lightweight and Easily Transported Consists of one infant manikin, size: 20" x 14" x 7", 10 lung/mouth protection bags, one insertion tool, and an instruction manual.

Q: What age patient does the Baby Buddy™ manikin simulate? A: It simulates a full-body infant, appropriate for neonatal and infant CPR training covering the two-finger chest compression technique,

correct hand placement landmarks, and rescue breathing.

Q: Why is a 5-pack more effective than a single manikin for group training? A: With one manikin per student, every trainee gets maximum hands-on practice time without waiting. Research consistently shows that compression quality and retention improves with individual practice repetitions. In a class of 20 students, 4 packs of 5 gives every student their own manikin simultaneously — transforming BLS training efficiency.

Q: Is this manikin hygienic for group training? A: Yes. The one-piece disposable lung/mouth protection bag is replaced between students — completely eliminating cross-contamination risk without any sanitisation of the manikin itself. The removable chest panel makes bag replacement quick and straightforward.

Q: Can this be used with AED training? A: Yes, with an important caveat — only AED training pads specifically designed for foam body manikins should be used. Standard electrode training pads may not be compatible with the foam body construction.

Q: How many training sessions will the included bags last? A: The pack includes 50 bags (10 per manikin). For a class of 20 students using 2 packs, this covers 5 complete group sessions before replacement bags need to be ordered.

Q: Are replacement lung/mouth protection bags available? A: Yes. Replacement disposable lung/mouth bags for the Baby Buddy™ manikin are available as a consumable — contact SEM Trainers to order.

Q: Is this suitable for NMC-compliant neonatal BLS skills labs? A: Yes. The Baby Buddy™ 5-Pack fulfils the infant CPR manikin requirements for NMC-compliant paediatric and neonatal BLS skills stations in medical college and nursing institute simulation labs.

Q: Can all 5 manikins be stored and transported together? A: Yes. All 5 manikins fit inside the included rugged carry bag, making the complete pack easy to move between training rooms, venues, or storage locations.

NOTE: Only use AED pads designed for foam body manikins.