



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Infant Choking Manikin

Item No. 1005778 [W44685]
Weight 2.7kg
MPN: 100-1640U / PP01640U / SB40035U

[Read More](#)

SKU:

Categories:BLS Child

Product Description

This 9-month-old baby brings to you the challenges associated with foreign body obstruction in the young pediatric age group. By placing the included obstruction in the pharyngeal cavity, you will occlude the throat and initialize the need for protocol maneuvers to remove the object. When proper hand placement occurs with sufficient force, the object expels from the mouth. Delivery with two objects, clothes and carry bag.

Item No. 1005778 [W44685], 2.7 kg, MPN 100-1640U / PP01640U / SB40035U, Nasco Healthcare, USA — Make field absent. 9-month-old infant. Same expelled-object consequence feedback mechanism as the Child (W44584) and Adolescent (W44583) choking torsos done earlier — completing the three-piece choking training family. Key clinical distinction: infant choking response

uses back blows and chest thrusts, not abdominal thrusts. The object is placed in the pharyngeal cavity (not behind the teeth as in older patients), making the mechanism and technique explicitly different. Supplied with two objects, clothes, and carry bag.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Infant Choking Manikin for Your Simulation Lab in India?

The Infant Choking Manikin (Item No. 1005778 / W44685, MPN: 100-1640U / PP01640U / SB40035U) by Nasco Healthcare, USA is a 2.7 kg 9-month-old infant choking response training manikin — the only patient age in the Nasco choking torso family that teaches the infant-specific protocol combining back blows and chest thrusts, which differs fundamentally from the abdominal thrust technique used for child and adolescent patients. The consequence-feedback mechanism is direct: an obstruction is placed in the pharyngeal cavity, and only when the trainee applies correct back blow or chest thrust technique with sufficient force is the object expelled — providing immediate, unambiguous pass/fail confirmation of correct performance. Supplied with two obstruction objects, clothes, and a carry bag. For schools, paediatric first-aid certification centres, NRP training programmes, nursing institutes, and community CPR outreach in India that need to teach infant choking response with built-in technique feedback at the correct 9-month-old patient proportions, this is the correct purpose-built platform.

The pharyngeal cavity placement of the obstruction — distinct from the subdiaphragmatic position that drives abdominal thrust training on child and adolescent torsos — accurately represents the clinical reality of infant airway obstruction where back blows are delivered with the infant face-down over the forearm, and chest thrusts are delivered with the infant face-up, using two fingers on the lower sternum rather than the encircling-hands abdominal grip used for older patients. The 2.7 kg weight and 9-month-old proportions correctly represent the held-infant position that back blow and chest thrust technique requires — where the trainee holds the manikin supported on the forearm and applies the blows or thrusts from behind. Compared to the Child Choking Torso (1005751/W44584, 2.08 kg) and Adolescent Choking Torso (1005750/W44583, 4.1 kg), this infant manikin presents the most technically distinct choking response protocol and is the correct companion for multi-age choking response curricula covering infant, child, and adolescent presentations.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: Does this manikin teach abdominal thrusts for infant choking? A: No. Infant choking response uses back blows and chest thrusts — not abdominal thrusts, which are contraindicated in infants. This manikin is designed to teach the correct infant-specific protocol: back blows with the infant face-down over the forearm and chest thrusts with the infant face-up.

Q: What is included with this manikin? A: Two obstruction objects, clothes, and a carry bag.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.