



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Life/form® Child CPR/Airway Management Torso

Item No. LF03633U

Weight 3.856kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Airway Management Child

Product Description

The CPR/Airway Management Torso is ideal for practicing intubation, ventilation, suctioning techniques, and CPR. This torso is designed to be the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest skin can be added later to provide IV access, intraosseous infusions, blood pressure monitoring, and ECG recognition and defibrillation training. Includes pump lubricant.

The Life/form® Child CPR/Airway Management Torso (LF03633U) is a paediatric simulation trainer designed for hands-on practice of CPR and airway management skills in a realistic 5-year-old child torso form. Built by Nasco Healthcare, USA, this trainer is the foundational unit of the Child CRiSis™ modular system — allowing institutions to start with core CPR and airway skills and expand later with

arms, legs, and defibrillation chest skin for full paediatric emergency simulation.

Q: What age group does this manikin simulate? A: It simulates a 5-year-old child, making it appropriate for paediatric CPR and airway management training scenarios.

Q: Can this torso be used with supraglottic airway devices? A: Yes. The manikin is fully compatible with supraglottic airway devices including LMA (Laryngeal Mask Airway) and bag-valve-mask (BVM) ventilation.

Q: Is this a standalone product or part of a larger system? A: It is both. It functions as a standalone CPR/airway trainer and also serves as the base unit of the Child CRiSis™ modular system. Arms, legs, and a defibrillation chest skin can be added later to enable IV access, intraosseous infusion, BP monitoring, and defibrillation training.

Q: What skills can be practised on this torso? A: CPR compressions, intubation, ventilation, suctioning, and supraglottic airway device insertion.

Q: Is this product suitable for NMC-compliant skills labs? A: Yes. It meets the simulation equipment standards required for paediatric emergency skills stations in NMC-compliant medical college skills labs.

Q: What is included in the package? A: The Child CPR/Airway Management Torso and pump lubricant.

Q: Who manufactures this product? A: It is manufactured by Nasco Healthcare, USA — a globally recognised brand in medical simulation equipment.

Q: Can I upgrade this unit to a full-body manikin later? A: Yes. The Child CRiSis™ system is modular. You can add limbs and additional chest skins over time without replacing the core torso unit.

Key Features:

- Realistic paediatric torso anatomy for accurate CPR compression training
- Compatible with supraglottic airway devices (LMA, BVM, etc.) for airway management practice
- Modular design — expandable to full Child CRiSis™ manikin with IV access, IO infusions, BP monitoring, and ECG/defibrillation training
- Supports intubation, ventilation, and suctioning technique practice
- Includes pump lubricant for maintenance
- Lightweight at 3.856 kg — easy to store and transport between lab
- **Ideal For:**
- Medical colleges and nursing institutes setting up NMC/INC-compliant paediatric skills labsPALS (Paediatric Advanced Life Support) training programmes

- Emergency medicine, paediatric nursing, and paramedic training
- AYUSH and allied health institutions with paediatric simulation requirements

Manikin can be used with Supraglottic airway

devices for airway training.

	Head/torso	Age	Digital intubation	Tongue edema	Laryngospasm	Stomach rise	Suctioning	Manual carotid pulse	Defibrillation	Additional features
NG Tube Insertion Trainer Item No: 1017250	head mounted on a board	1-year old	-	-	-	Yes	-	-	-	trainer for instruction of nasogastric tube insertion and care
Infant Airway Simulator Item No: 1018863	full-body	1-year old	-	-	-	-	Yes	-	-	chest cavity containing realistic organs, nasal passage permits placement of NP tube
Infant Airway Trainer Item No: 1017954	head mounted on a board	1-year old	-	-	-	Yes	Yes	-	-	accepts 4.0 mm endotracheal tubes
Infant Intubation Head Item No: 1017245	head mounted on a board	1-year old	-	-	-	-	Yes	-	-	accepts 2.5 mm endotracheal tubes
Advanced Infant Intubation Head Item No: 1017236	head mounted on a board	1-year old	-	Yes	-	Yes	Yes	-	-	Vocal cords are highlighted for easy viewing
Airway Management Trainer Item No: 1014201	head mounted on a board	3-year old	-	-	-	-	Yes	-	-	accepts 4.0 mm endotracheal tubes

Advanced Child Airway Management Trainer <i>Item No:</i> 1019787	head	5-year old	-	Yes	Yes	-	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
Child Airway Management Trainer Head Item <i>No:</i> 1019788	head	5-year old	-	-	-	Yes	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
Advanced Child Airway Management Trainer <i>Item No:</i> 1005698	head mounted on a board	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	-	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
Basic Child CPR/ Airway Management Torso" <i>Item No:</i> 1018865	torso	5-year old	-	-	-	-	-	-	-	This torso is the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest can be added later.

Advanced Child CPR/Airway Management Torso <i>Item No:</i> 1018864	torso	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	Yes	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
---	-------	------------	-----	-----	-----	-----	-----	-----	-----	---