



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Infant Airway Trainer with Stand

Item No. LF03623U
Weight 6.4kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Airway Management Child



Product Description

The infant airway management trainer features realistic size and anatomy and is perfect for practicing airway management skills on an infant. Realistic anatomy includes mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea, and esophagus. Separate lungs for auscultation. Inflatable stomach bladder indicated esophageal insufflation. Practice suctioning and oral intubation techniques. Each

trainer includes a stand to mount head, pump spray lubricant, and hard carrying case. Can be attached to Basic Infant CRiSis™ manikin (W44717/1017248)

Why Choose the Life/form® Infant Airway Trainer with Stand for Your Simulation Lab in India?

The Life/form® Infant Airway Trainer with Stand (Item No. LF03623U / 1017954) by Nasco Healthcare, USA is a 6.4 kg infant airway management training head mounted on a stand — the complete standalone infant airway training platform in the Life/form® infant airway range, with a hard carry case and pump spray lubricant included as standard. A downloadable user manual is available on the product page. From the comparison table, this configuration supports stomach rise and suctioning, and accepts 4.0 mm endotracheal tubes — the correct paediatric ET tube size for a 1-year-old patient simulation. For NRP training programmes, nursing institutes, paediatric emergency medicine skills labs, and NMC-mandated simulation labs in India that need a standalone infant airway management platform with stand and carry case for regular training and outreach sessions, this is the correct ready-to-use configuration.

The anatomical landmark set covers mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea, and oesophagus — the complete reference set for oral intubation technique training at 1-year-old infant proportions. Separate bifurcated lungs enable bilateral auscultation confirmation of correct tube placement. The inflatable stomach bladder provides immediate oesophageal intubation consequence feedback — stomach rises instead of lungs when the tube is incorrectly placed. Suctioning technique training is supported. The stand presents the infant head at the correct approach angle for infant intubation technique development. The head can also be detached from the stand and attached directly to the Basic Infant CRiSis™ Manikin (W44717 / 1017248) for combined infant CPR and airway management training on the same platform. The hard carry case makes this the most transport-ready infant airway trainer in the Nasco range, suitable for outreach training and multi-venue programmes.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: What ET tube size does this trainer accept? A: The Life/form® Infant Airway Trainer with Stand accepts 4.0 mm endotracheal tubes — the appropriate size for a 1-year-old infant patient simulation.

Q: Can this head be used with the Basic Infant CRiSis™ Manikin? A: Yes. The head can be detached from the stand and attached directly to the Basic Infant CRiSis™ Manikin (W44717 /

1017248) for combined CPR and airway training.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. A downloadable user manual is available on the product page. For support, reach us at sem@semtrainers.com or call +91-88495 63724.

	Head/torso	Age	Digital intubation	Tongue edema	Laryngospasm	Stomach rise	Suctioning	Manual carotid pulse	Defibrillation	Additional features
NG Tube Insertion Trainer Item No: 1017250	head mounted on a board	1-year old	-	-	-	Yes	-	-	-	trainer for instruction of nasogastric tube insertion and care
Infant Airway Simulator Item No: 1018863	full-body	1-year old	-	-	-	-	Yes	-	-	chest cavity containing realistic organs, nasal passage permits placement of NP tube
Infant Airway Trainer Item No: 1017954	head mounted on a board	1-year old	-	-	-	Yes	Yes	-	-	accepts 4.0 mm endotracheal tubes
Infant Intubation Head Item No: 1017245	head mounted on a board	1-year old	-	-	-	-	Yes	-	-	accepts 2.5 mm endotracheal tubes
Advanced Infant Intubation Head Item No: 1017236	head mounted on a board	1-year old	-	Yes	-	Yes	Yes	-	-	Vocal cords are highlighted for easy viewing
Airway Management Trainer Item No: 1014201	head mounted on a board	3-year old	-	-	-	-	Yes	-	-	accepts 4.0 mm endotracheal tubes

Advanced Child Airway Management Trainer <i>Item No:</i> 1019787	head	5-year old	-	Yes	Yes	-	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
Child Airway Management Trainer Head Item <i>No:</i> 1019788	head	5-year old	-	-	-	Yes	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
Advanced Child Airway Management Trainer <i>Item No:</i> 1005698	head mounted on a board	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	-	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
Basic Child CPR/ Airway Management Torso" <i>Item No:</i> 1018865	torso	5-year old	-	-	-	-	-	-	-	This torso is the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest can be added later.

Advanced Child CPR/Airway Management Torso <i>Item No:</i> 1018864	torso	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	Yes	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
---	-------	------------	-----	-----	-----	-----	-----	-----	-----	---