




## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Advanced Infant Intubation Head with Board

Item No. 101-130U  
Weight 4.08kg  
Dimensions 38.1 x 25.4 x 20.3 cm  
Make: Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:**Airway Management Child

## Product Description

This station trainer incorporates new skin technology. The advantages of durability and lifelike appearance, in association with the convenience of a lightweight stand, put this device at the head of the class for pediatric airway training. This new material eliminates tearing of the airway, saving the need for costly repairs caused by beginning student intubation attempts. Translucent property allows for lifelike illumination of the airway and neck as the skill is attempted. Vocal cords are highlighted for easy viewing, the tongue swells, and all the anatomical landmarks are present.

**Why Choose the Life/form® Advanced Child CPR/Airway Management Torso with Defibrillation Features for Your Simulation Lab in India?**

The Life/form® Advanced Child CPR/Airway Management Torso with Defibrillation Features (Item No. LF03763U / 1018864) by Nasco Healthcare, USA is an 8.1 kg paediatric torso simulating a 5-year-old

patient — and the only configuration in the Life/form® Child Airway Management Trainer range to combine all seven advanced training capabilities on a single platform: digital intubation, tongue oedema, laryngospasm, stomach rise, suctioning, manual carotid pulse, and defibrillation. It is also the foundation torso of the Child CRiSis™ system — arms and legs for IV access, intraosseous infusion, and blood pressure monitoring can be added at any time, making this the correct choice for PALS simulation labs in India that want a scalable platform that grows with their training programme. For simulation labs that need to run paediatric airway management, difficult airway scenarios, CPR, and ACLS defibrillation training on a single torso without a separate defibrillation trainer, this is purpose-built for that requirement.

The defibrillation chest skin incorporates an internal load box that absorbs the full strength of every shock from any standard manual, automatic, or semi-automatic external monitor defibrillator — protecting both trainees and equipment. Four ECG monitoring sites and two defibrillation pad sites are present, enabling ECG monitoring and shock delivery training with the trainee's own department defibrillator. The Optional Interactive ECG Simulator (LF03670U) can be added to enable ECG rhythm recognition and arrhythmia training through the same torso. The airway module — with slightly anterior position, instructor-activated tongue swelling and laryngospasm — supports digital, oral, and nasal intubation alongside ET, EOA, PTL, LMA, EGTA, Combitube®, King Systems, and all supraglottic airway devices. Jaw thrust technique and Sellick manoeuvre training are supported. Separate bifurcated lungs enable bilateral auscultation confirmation of correct tube placement. The inflatable stomach bladder provides oesophageal intubation consequence feedback. CPR training uses the standard paediatric hand placement and compression depth for a 5-year-old torso. Manual carotid pulse enables integrated assessment within the paediatric resuscitation sequence.

**Why Buy from SEM Trainers?** SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

## Frequently Asked Questions

**Q: Can arms and legs be added to this torso later?** A: Yes. LF03763U is the foundation torso of the Life/form® Child CRiSis™ system. IV access arms, IO infusion legs, and blood pressure monitoring arms can be added at any time to expand procedural training capability without replacing the torso.

**Q: Is an ECG Simulator included?** A: No. The defibrillation chest supports ECG monitoring and shock delivery with any standard defibrillator brand. For ECG rhythm recognition and arrhythmia training, the optional Interactive ECG Simulator (LF03670U) can be added separately.

**Q: What is the warranty and how do I get support?** A: All Nasco Healthcare products carry standard

manufacturer warranty. For support, reach us at [sem@semtrainers.com](mailto:sem@semtrainers.com) or call +91-88495 63724.

|                                                                            | Head/torso              | Age        | Digital intubation | Tongue edema | Laryngospasm | Stomach rise | Suctioning | Manual carotid pulse | Defibrillation | Additional features                                                                  |
|----------------------------------------------------------------------------|-------------------------|------------|--------------------|--------------|--------------|--------------|------------|----------------------|----------------|--------------------------------------------------------------------------------------|
| <b>NG Tube Insertion Trainer</b><br>Item No: <a href="#">1017250</a>       | head mounted on a board | 1-year old | -                  | -            | -            | Yes          | -          | -                    | -              | trainer for instruction of nasogastric tube insertion and care                       |
| <b>Infant Airway Simulator</b><br>Item No: <a href="#">1018863</a>         | full-body               | 1-year old | -                  | -            | -            | -            | Yes        | -                    | -              | chest cavity containing realistic organs, nasal passage permits placement of NP tube |
| <b>Infant Airway Trainer</b><br>Item No: <a href="#">1017954</a>           | head mounted on a board | 1-year old | -                  | -            | -            | Yes          | Yes        | -                    | -              | accepts 4.0 mm endotracheal tubes                                                    |
| <b>Infant Intubation Head</b><br>Item No: <a href="#">1017245</a>          | head mounted on a board | 1-year old | -                  | -            | -            | -            | Yes        | -                    | -              | accepts 2.5 mm endotracheal tubes                                                    |
| <b>Advanced Infant Intubation Head</b><br>Item No: <a href="#">1017236</a> | head mounted on a board | 1-year old | -                  | Yes          | -            | Yes          | Yes        | -                    | -              | Vocal cords are highlighted for easy viewing                                         |
| <b>Airway Management Trainer</b><br>Item No: <a href="#">1014201</a>       | head mounted on a board | 3-year old | -                  | -            | -            | -            | Yes        | -                    | -              | accepts 4.0 mm endotracheal tubes                                                    |

|                                                                                                |                         |            |     |     |     |     |     |     |   |                                                                                                                        |
|------------------------------------------------------------------------------------------------|-------------------------|------------|-----|-----|-----|-----|-----|-----|---|------------------------------------------------------------------------------------------------------------------------|
| <b>Advanced Child Airway Management Trainer</b><br><i>Item No:</i><br><a href="#">1019787</a>  | head                    | 5-year old | -   | Yes | Yes | -   | Yes | Yes | - | can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins                     |
| <b>Child Airway Management Trainer Head Item</b><br><i>No:</i> <a href="#">1019788</a>         | head                    | 5-year old | -   | -   | -   | Yes | Yes | Yes | - | can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins                     |
| <b>Advanced Child Airway Management Trainer</b><br><i>Item No:</i><br><a href="#">1005698</a>  | head mounted on a board | 5-year old | Yes | Yes | Yes | Yes | Yes | Yes | - | also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion                          |
| <b>Basic Child CPR/ Airway Management Torso"</b><br><i>Item No:</i><br><a href="#">1018865</a> | torso                   | 5-year old | -   | -   | -   | -   | -   | -   | - | This torso is the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest can be added later. |

|                                                                                                 |       |            |     |     |     |     |     |     |     |                                                                                               |
|-------------------------------------------------------------------------------------------------|-------|------------|-----|-----|-----|-----|-----|-----|-----|-----------------------------------------------------------------------------------------------|
| <b>Advanced Child CPR/Airway Management Torso</b><br><i>Item No:</i><br><a href="#">1018864</a> | torso | 5-year old | Yes | Yes | Yes | Yes | Yes | Yes | Yes | also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion |
|-------------------------------------------------------------------------------------------------|-------|------------|-----|-----|-----|-----|-----|-----|-----|-----------------------------------------------------------------------------------------------|