



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

AirSim Child Combo Bronchi X

Item No. 1021903
Weight 2.612 kg
MPN: CC50006X

[Read More](#)

SKU:

Categories: Airway Management Child



Product Description

The AirSim Child Combo Bronchi X is the 'all in-one' manikin for pediatric airway management. It is ideal for practising oral and nasal intubation training techniques, percutaneous tracheostomy, emergency cricothyroidotomy and bronchoscopy techniques on a six year old child.

Features include:

- Innovative pediatric AirSim® airway and nasal cavity created from DT DICOM data
- Anatomically correct internal features and visually accurate landmarks in airway and nasal passage
- Wraparound replaceable neck skin rotates for up to 20 incisions
- Easy identification of sternal notch, tracheal rings and cricoid and laryngeal cartilages
- Inflatable tongue with lifelike texture can simulate edema
- Silicone skin covering with lifelike texture
- Replaceable lung bags inflate with successful ventilation
- Optional difficult airways available for more advanced training
- Optional subcutaneous fat tissue demonstrates more difficult tracheostomy and cricothyroidotomy procedures

The AirSim Child Combo Bronchi X will provide a solution to your training needs in:

- Needle and surgical cricothyroidotomy
- Percutaneous tracheostomy
- Practising oral and nasal intubation techniques
- Practising laryngeal mask airway and supraglottic device insertion
- Practising combi tube insertion
- Practising double lumen endotracheal tube insertion
- Practising fibre optic intubation and video laryngoscopy
- Practising bag mask ventilation techniques
- Diagnostic bronchoscopy
- Practising lung isolation techniques using left and right endobronchial tubes and bronchial blockers
- Practising lung suctioning techniques

Product components:

- 1 Carrier/storage case
- 1 USB pen drive user manual
- 1 Bottle of TruCorp® lubrication
- 3 Adult lung bags
- 5 Larynx inserts (1 attached, 4 extra)
- 1 Wrap around neck skin (attached)

Manikin can be used with Supraglottic airway devices for airway training.