



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Advanced KERi™ Nursing Manikin

Item No. 1017567 [W44403]

Weight 12.6 kg

Dimensions 94 x 53.4 x 27.9 cm

MPN: LF04022U

[Read More](#)

SKU:

Categories:Adult Patient Care

Product Description

KERi™ has a non-age-specific appearance. Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - moves like a person.

- Lightweight - approximately 28 lbs.
- Full-size adult manikin - measures 58"
- Overall female appearance
- Converts to male with removal of wig and attachment of male genitals

Visual Inspection:

- Normal and Cancerous Mole
- Stage 1 Sacral Ulcer
- Dilated and constricted pupils
- Reddened Skin Folds

Patient Care Simulation:

- Bandaging and Wound Dressing
- Bed Baths
- Clothing Changes
- Denture Placement and Removal, Upper and Lower
- Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement
- Eye Irrigation
- Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded
- Hair Care Washing and Combing
- Intramuscular Injection Sites: Arm, Thigh, Buttock
- Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning
- Ostomy Care: Ileostomy and Colostomy Tissue
- Maintenance and Appliance Application, Lavage, and Suctioning
- Patient Positioning
- Patient Transfer Techniques
- Pericare
- Tracheostomy Care, Lavage and Suctioning
- Gastrostomy Procedures, Lavage and Gavage
- Enema Administration, Female Only
- Urinary Catheterization, Female and Uncircumcised Male
- Pap Smears and Douching
- Prostate Exam, Stage B
- Nasogastric Tube Placement
- Internal Fluid Reservoirs

The Advanced KERi™ Manikin simulates over 35 nursing and medical procedures, and also includes:

- Right IV Training Arm with a shoulder intramuscular injection site featuring replaceable skin and veins that roll when palpated; realistic flashback confirms proper needle placement
- Left Blood Pressure Training Arm that reproduces all five Korotkoff sounds and allows the instructor to vary systolic and diastolic levels, pulse rate, volume, and auscultatory gap

- Standard arms with right and left shoulder intramuscular injection sites are also included.

	Basic: KERi™ 1020187 GERi™ 1005620	Complete: KERi™ 1013742 GERi™ 1005597	Advanced: KERi™ 1017567 GERi™ 1005606	Auscultation: KERi™ 1020147 GERi™ 1020146
Skin inspection and care	Yes	Yes	Yes	Yes
Bandaging	Yes	Yes	Yes	Yes
Washing	Yes	Yes	Yes	Yes
Denture Placement	Yes	Yes	Yes	Yes
Hearing aid placement	Yes	Yes	Yes	Yes
Ear and eye irrigation	Yes	Yes	Yes	Yes
Ostomy and tracheostomy care	Yes	Yes	Yes	Yes
Positioning and transfer techniques	Yes	Yes	Yes	Yes
Pericare	Yes	Yes	Yes	Yes
Intramuscular injections	Yes	Yes	Yes	Yes
Gastrostomy	-	Yes	Yes	Yes
Pap smears and douching	-	Yes	Yes	Yes
Enema administration	-	Yes	Yes	Yes
Prostate exam	-	Yes	Yes	Yes
Urinary catheterisation	-	Yes	Yes	Yes
Nasogastric tube placement	-	Yes	Yes	Yes
IV insertion	-	-	Yes	
Blood pressure monitoring	-	-	Yes	
Recognition of heart sounds	-	-	-	Yes
Recognition of lung sounds	-	-	-	Yes