



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

GERi™ Auscultation Manikin

Item No. 1020146
Weight 14.8kg
MPN: LF04003U

[Read More](#)

SKU:

Categories: Geriatric Patient Care



Product Description

GERi™ has an elderly appearance with skin wrinkles and folds. Complete lifelike range of motion,

realistic patient positioning, and non-pinching joints - moves like a person.

- Lightweight - approximately 28 lbs. / ca. 12,7 kg
- Full-size adult manikin - measures 58" / 1,47 m
- Overall female appearance
- Converts to female with removal of wig and attachment of female genitals

Visual Inspection:

- Normal and Cancerous Mole
- Stage 1 Sacral Ulcer
- Dilated and Constricted Pupils
- Reddened Skin Folds

Patient Care Simulation:

- Bandaging and Wound Dressing
- Bed Baths
- Clothing Changes
- Denture Placement and Removal, Upper and Lower
- Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement
- Eye Irrigation
- Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded
- Hair Care Washing and Combing
- Intramuscular Injection Sites: Arm, Thigh, Buttock
- Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning
- Ostomy Care: Ileostomy and Colostomy Tissue
- Maintenance and Appliance Application, Lavage, and Suctioning
- Patient Positioning
- Patient Transfer Techniques
- Pericare
- Tracheostomy Care, Lavage and Suctioning
- Gastrostomy Procedures, Lavage and Gavage
- Enema Administration, Female Only
- Urinary Catheterization, Female and Uncircumcised Male
- Pap Smears and Douching
- Prostate Exam, Stage B
- Nasogastric Tube Placement
- Internal Fluid Reservoirs

The GERi™ Auscultation Manikin also includes:

- Auscultation Sites with Heart and Lung Sounds

- Heart Sites, Six anterior locations with 12 heart conditions
- Heart Conditions: • Aortic Regurgitation, Atrial Septal Defect, Holosystolic, Midsystolic, Mitral Stenosis, Normal, PDA, Pulmonary Stenosis, S3 Gallop, S4 Gallop, Systolic Click, and VSD
- Lung Sites: Five Anterior, six upper posterior, four lower posterior, and two midaxillary locations with 12 lung conditions
- Lung Conditions: Bronchial, Cavernous, Coarse Crackle, Fine Crackle, Friction Rub, Mono Wheeze, Normal, Pectoriloquy, Pulmonary Edema, Rhonchi, Stridor, and Wheeze
- **Simulates over 35 nursing and medical procedures**

	Basic: KERi™ 1020187 GERi™ 1005620	Complete: KERi™ 1013742 GERi™ 1005597	Advanced: KERi™ 1017567 GERi™ 1005606	Auscultation: KERi™ 1020147 GERi™ 1020146
Skin inspection and care	Yes	Yes	Yes	Yes
Bandaging	Yes	Yes	Yes	Yes
Washing	Yes	Yes	Yes	Yes
Denture Placement	Yes	Yes	Yes	Yes
Hearing aid placement	Yes	Yes	Yes	Yes
Ear and eye irrigation	Yes	Yes	Yes	Yes
Ostomy and tracheostomy care	Yes	Yes	Yes	Yes
Positioning and transfer techniques	Yes	Yes	Yes	Yes
Pericare	Yes	Yes	Yes	Yes
Intramuscular injections	Yes	Yes	Yes	Yes
Gastrostomy	-	Yes	Yes	Yes
Pap smears and douching	-	Yes	Yes	Yes
Enema administration	-	Yes	Yes	Yes
Prostate exam	-	Yes	Yes	Yes
Urinary catheterisation	-	Yes	Yes	Yes
Nasogastric tube placement	-	Yes	Yes	Yes
IV insertion	-	-	Yes	
Blood pressure monitoring	-	-	Yes	
Recognition of heart sounds	-	-	-	Yes
Recognition of lung sounds	-	-	-	Yes

GERi™ Auscultation Manikin (1020146 / LF04003U)

GERi™ has an elderly appearance with skin wrinkles and folds, offering complete lifelike range of motion, realistic patient positioning, and non-pinching joints that move like a person.

A true full-size, lightweight adult manikin. GERi™ measures 58" and weighs approximately 28 lbs. (~12.7 kg), with an overall female appearance that converts to male with removal of the wig and attachment of male genitals.

Visual inspection findings built right in: normal and cancerous mole, a stage 1 sacral ulcer, dilated and constricted pupils, and reddened skin folds.

A comprehensive range of patient care simulation skills, including bandaging and wound dressing, bed baths, clothing changes, denture placement and removal, ear canal irrigation with otic drops and hearing aid placement, eye irrigation, finger and toe manipulation with bandaging, hair care, intramuscular injection sites (arm, thigh, buttock), oral and nasal hygiene (lavage, gavage, suctioning), ostomy care with appliance maintenance, patient positioning and transfer techniques, pericare, tracheostomy care, gastrostomy procedures, enema administration (female), urinary catheterization (female and uncircumcised male), pap smears and douching, prostate exam (Stage B), nasogastric tube placement, and internal fluid reservoirs.

Auscultation tier adds comprehensive heart and lung sound training. Six anterior heart sites cover 12 heart conditions: aortic regurgitation, atrial septal defect, holosystolic, midsystolic, mitral stenosis, normal, PDA, pulmonary stenosis, S3 gallop, S4 gallop, systolic click, and VSD. Lung sites — five anterior, six upper posterior, four lower posterior, and two midaxillary locations — cover 12 lung conditions: bronchial, cavernous, coarse crackle, fine crackle, friction rub, mono wheeze, normal, pectoriloquy, pulmonary edema, rhonchi, stridor, and wheeze.

Simulates over 35 nursing and medical procedures.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical education resources to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: How many heart and lung conditions can be auscultated? A: 12 heart conditions across 6 anterior sites, and 12 lung conditions across 17 total lung sites.

Q: Can this manikin be configured as male? A: Yes — the overall appearance is female, but it converts

to male with removal of the wig and attachment of male genitals.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.