




SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Vein Manifold for 1005678

Item No. LF00993U
Weight 0.078kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:BLS and CPR Accessories

Product Description

Replacement for use with Life/form® Advanced Venipuncture and Injection Arms. The Vein Manifold for 1005678 (Item No. LF00993U) is a genuine replacement component manufactured by Nasco Healthcare, USA, designed specifically for use with the Life/form® Advanced Venipuncture and Injection Arms. Weighing just 0.078 kg, this lightweight replacement vein manifold is an essential spare part for medical training institutions, nursing institutes, and clinical skills labs that regularly use Life/form® Advanced Venipuncture and Injection Arm trainers for IV cannulation and blood draw training. Using a genuine replacement vein manifold ensures continued proper functionality, realistic blood flow simulation, and accurate training across multiple sessions.

Q1. What is the Vein Manifold for 1005678 used for? It is a replacement vein manifold component designed for use with the Life/form® Advanced Venipuncture and Injection Arms, helping maintain

realistic and accurate blood flow simulation during IV cannulation and venipuncture training sessions.

Q2. Which trainer is this vein manifold compatible with? This vein manifold is designed as a genuine replacement for the Life/form® Advanced Venipuncture and Injection Arm (model 1005678). Please verify compatibility with your specific trainer model before ordering.

Q3. What is a vein manifold and why does it need replacing? The vein manifold is an internal component that distributes artificial blood flow to the multiple vein sites on the injection arm trainer. After repeated use, it can wear out or develop leaks, reducing the realism and effectiveness of venipuncture training — making timely replacement essential.

Q4. Why should I use a genuine replacement vein manifold? Using a genuine Life/form® replacement component ensures a proper fit and reliable performance, maintaining the accuracy and longevity of the Advanced Venipuncture and Injection Arm without the risk of damage from incompatible third-party parts.

Q5. What is the weight of this product? The vein manifold weighs just 0.078 kg, making it compact and easy to store alongside the training arm.

Q6. Who manufactures this product? It is manufactured by Nasco Healthcare, USA — a globally trusted brand in medical simulation and clinical skills training equipment.

Q7. Who is this product ideal for? It is ideal for medical colleges, nursing institutes, clinical skills labs, and healthcare training centers that regularly use the Life/form® Advanced Venipuncture and Injection Arm for IV access and blood draw training.

Q8. Where can I purchase this in India? You can enquire and order through SEM Trainers at **+91-88495 63724** or **+91-98791 03905**, or visit semtrainers.com.