



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Arm/hand vein set for IV injection hand

Item No. 140-143

Weight 1.543kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

Arm/hand vein set for IV injection hand W44600

Life/form® Arm/Hand Vein Set for IV Injection Hand W44600 (140-143) — Genuine OEM Replacement Consumable

The Arm/Hand Vein Set for IV Injection Hand (Item No. 140-143) is a genuine Nasco Healthcare replacement vein assembly for the **Life/form® IV Injection Hand trainer (W44600)**. It is the correct reorder item when the vein tubing in your W44600 trainer has been exhausted through repeated IV cannulation and venipuncture practice sessions — restoring the trainer's simulated blood flashback function and vein palpation fidelity without replacing the entire unit.

The replacement vein set includes the complete internal vein tubing assembly for the W44600 arm and hand trainer. Once installed, the trainer delivers consistent needle entry feedback, palpable vein resistance, and simulated blood flashback confirmation matching the original factory specification. Using a genuine Nasco Healthcare OEM vein set ensures correct tubing diameter, wall thickness, and connector fit for the W44600 platform.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical simulators and genuine spare parts to 1,500+ institutions across 28 states in India. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which trainer is this vein set compatible with? A: This vein set is the OEM replacement for the Life/form® IV Injection Hand trainer — Nasco Healthcare reference W44600 only. It is not compatible with other Nasco Healthcare IV arm or hand trainers, including the Advanced IV Hand (W44117) or the Advanced Venipuncture and Injection Arm (LF01121U / W44216).

Q: Is this a genuine OEM part? A: Yes. Item No. 140-143 is a genuine Nasco Healthcare Life/form® replacement consumable.

Q: What is the difference between this part and the Advanced IV Hand Replacement Skin and Veins (LF01140U)? A: These are replacement parts for two different parent trainers. LF01140U is the skin and vein replacement for the Advanced IV Hand (W44117). This item (140-143) is the vein set replacement for the IV Injection Hand (W44600). Confirm your trainer's model reference before ordering.

Q: When should the vein set be replaced? A: When the vein tubing no longer produces consistent simulated blood flashback on correct needle placement, or when the tubing wall has been perforated to the point that it no longer holds fluid pressure reliably between sessions.