



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Liquid Skin Repair for Intradermal Injection Simulator

Item No. LF01009U

Weight 0.088kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

2-oz. bottle. The liquid skin repair will extend the life of skin. For use on **Life/form®** replicas ONLY. **Not for human use.**

Liquid Skin Repair for Life/form® Intradermal Injection Simulator (LF01009U)

The Liquid Skin Repair (Item No. LF01009U) is the official Nasco Healthcare / Life/form® consumable maintenance fluid for the Life/form® Intradermal Injection Simulator (W44097 / LF01008U) — the forearm ID injection simulator whose skin forms a visible welt after each intradermal injection, replicating the clinical confirmation of correct ID placement. Over time and with repeated use, these welt marks accumulate and the skin surface requires restoration. The Liquid Skin Repair fills and seals existing needle marks and welt formations, extending the usable training life of the simulator skin

without requiring a full skin replacement.

Compatible with Life/form® Intradermal Injection Simulator skin only. Not compatible with other Life/form® or Nasco Healthcare simulators. **Not for human use.**

Supplied as: one 2 oz (59 ml) bottle — sufficient for multiple skin restoration treatments depending on usage volume.

Order alongside the Intradermal Injection Simulator (W44097) or as a standalone consumable replenishment. Contact SEM Trainers for pricing and availability of the Replacement Skin Set for the W44097 simulator (available separately when liquid repair is no longer sufficient).

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical simulators and consumables to 1,500+ institutions across 28 states in India. We provide pan-India delivery and dedicated after-sales support.

Q: Which simulator is this consumable for? A: The Life/form® Intradermal Injection Simulator (Item No. W44097 / LF01008U) only. It is not compatible with other Life/form® or Nasco Healthcare simulators.

Q: What does the liquid skin repair do? A: It fills and seals welt marks and needle punctures that accumulate on the simulator skin after repeated intradermal injection exercises, restoring the surface for continued training use and extending the life of the skin before a full replacement is needed.

Q: How much liquid is in the bottle? A: 2 oz (59 ml) per bottle.

Q: When should I order a replacement skin instead of liquid repair? A: When the skin surface is too damaged for liquid repair to restore adequately. Contact SEM Trainers to enquire about the Replacement Skin Set for the W44097 simulator.

Q: Is this safe for use on human skin? A: No. This product is for Life/form® simulator skin only. Not for human use.