



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Life/form® Infant IV Leg

Item No. LF03636U
Weight 0.5kg
Dimensions 0.001 x 0.001 x 0.001 cm
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

A special, extremely thin, synthetic skin is paired with rubber tubing featuring a small lumen and thin walls. The greater and lesser saphenous veins are accessible, as well as the dorsal venous arch on the foot. Comes with two IV bags with clamps, one pint Life/form® blood, one 3 cc syringe, one 22-gauge needle, and winged infusion set. Fluid stand not included.

Why Choose the Life/form® Infant IV Leg for Your Nursing Skills Lab in India?

The Life/form® Infant IV Leg (Item No. LF03636U, Weight: 0.5kg, Make: Nasco Healthcare, USA) gives students a dedicated way to practice infant IV access at the lower extremity — a standard alternative site when arm veins aren't accessible in real infant patients. Like its companion arm model, it uses an extremely thin synthetic skin paired with rubber tubing featuring a small lumen and thin walls, closely replicating the scale and fragility of real infant vasculature.

The greater and lesser saphenous veins are accessible, along with the dorsal venous arch on the foot — the key sites used for infant IV starts in the leg and foot.

What's included: Two IV bags with clamps, one pint of Life/form® simulated blood, one 3cc syringe, one 22-gauge needle, and a winged infusion set. (*Fluid stand not included.*)

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which veins can be accessed on this leg? A: The greater and lesser saphenous veins, along with the dorsal venous arch on the foot.

Q: Is this different from the Life/form® Infant IV Arm? A: Yes — this trainer focuses on lower-extremity access (leg and foot), while the Infant IV Arm covers upper-extremity sites. Both share the same thin-skin, small-lumen design for infant-scale realism.

Q: Is a fluid stand included? A: No — the fluid stand is sold separately.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.