



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



12 Intraosseous Bones with Blood Capsules

Item No. 101-092

Weight 0.5kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

Ship weight 2 lbs.

Life/form® 12 Intraosseous Bones with Blood Capsules (101-092) — Replacement Consumable Pack

The 12 Intraosseous Bones with Blood Capsules (Item No. 101-092) is a genuine Nasco Healthcare replacement consumable pack for Life/form® intraosseous infusion training simulators — providing 12 replacement IO bone inserts, each pre-loaded with a blood capsule that ruptures on correct needle placement to simulate the flashback confirmation used in live IO access procedures.

This pack is the correct reorder item when the IO bone inserts in your Life/form® IO simulator have

been exhausted through repeated practice sessions. Each bone insert is designed for single-use or limited-use training: the blood capsule confirms needle placement once, after which the insert should be replaced to maintain training fidelity for the next trainee. Ordering this pack restores your simulator to full training capacity without replacing the entire unit.

Compatible with Life/form® intraosseous infusion simulators in the Nasco Healthcare range, including the Life/form® Adult Sternal Intraosseous Infusion Simulator (available separately on this site). Confirm compatibility with your specific simulator model before ordering.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical simulators and genuine spare parts to 1,500+ institutions across 28 states in India. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What is included in this pack? A: 12 intraosseous bone inserts, each pre-loaded with a blood capsule for needle placement confirmation feedback.

Q: How does the blood capsule work? A: When the IO needle is correctly placed through the cortex of the bone insert, it ruptures the blood capsule inside — producing simulated blood flashback that confirms correct placement, exactly as the marrow blood return confirms correct placement in a live IO procedure.

Q: Are these inserts single-use? A: Yes — each bone insert and blood capsule is designed for single-use training to maintain consistent feedback fidelity. Once the blood capsule has been ruptured, the insert should be replaced.

Q: Which simulators are these compatible with? A: These replacement bones are for Life/form® intraosseous infusion simulators from Nasco Healthcare. Confirm compatibility with your specific model before ordering — contact SEM Trainers if unsure.

Q: How many training sessions does one pack support? A: One pack of 12 bones supports 12 individual IO needle placement practice attempts, or fewer if multiple attempts are made per trainee per session.