




## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Intradermal Injection Simulator

**Item No.** LF01008U  
**Weight** 0.77kg  
**Dimensions** 33 x 28 x 13 cm  
**Make:** Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:** Injections and Punctures

## Product Description

The simulator features a forearm from the wrist to just below the elbow. Vinyl skin provides realistic feel and appearance to ensure a realistic training experience. The simulator features eight sites for practicing intradermal injections. If fluid is properly injected, a characteristic skin welt will form. The welt is removed by withdrawing the fluid after practice. Each site is reusable by dozens of students. Supplied with sealant, syringe and storage box.

	Intramuscular Injection Simulator - Upper Arm	Intramuscular Injection Simulator - Upper Leg	Buttock Injection Simulator	Intramuscular injection simulator	Intramuscular Injection Model of Buttocks	Intramuscular Injection Simulator	Intradermal Injection Simulator
<b>Item No.</b>	<a href="#">1009840</a>	<a href="#">1000511</a>	<a href="#">1000514</a>	<a href="#">1010008</a>	<a href="#">1005394</a>	<a href="#">1005586</a>	<a href="#">1005628</a>
<b>Body part</b>	upper arm	upper leg	buttock	muscle extract with skin	pelvis with buttocks	lower torso	forearm

<b>Palpable anatomical landmarks</b>	✓	✓	✓	-	✓	✓	-
<b>Injection type</b>	intramuscular	intramuscular	intramuscular	intramuscular	intramuscular	intramuscular	intra dermal
<b>Injection sites</b>	as indicated by a trainer	as indicated by a trainer	as indicated by a trainer	all the surface	as indicated by a trainer	as indicated by a trainer	eight
<b>Proper injection</b>	audiovisual feedback	audiovisual feedback	audiovisual feedback	-	green light will appear	-	skin welt will form
<b>Additional information</b>	training and test modes are available	training and test modes are available	training and test modes are available	will absorb fluid	special skin material does not show needle marks	simulated bony structure is embedded in the torso	liquid skin repair kit included (M-1013686)

**Why Choose the Intradermal Injection Simulator for Your Simulation Lab in India?** The Intradermal Injection Simulator from 3B Scientific is a trusted training tool used across simulation labs in India for hands-on intradermal injection practice. This simulator helps medical students and healthcare professionals build critical competencies in a safe, controlled environment — fully aligned with NMC's Competency-Based Medical Education (CBME) guidelines for mandatory clinical skills training. **Who Is This Simulator Ideal For?** This product is recommended for medical colleges setting up OSCE stations, nursing institutions, and hospitals building clinical skills laboratories as per NMC guidelines. Whether you are setting up a new simulation lab or upgrading your existing facility, the Intradermal Injection Simulator offers an excellent balance of realism, durability, and training value. **Why Buy from SEM Trainers?** SEM Trainers is a leading supplier of medical manikins and simulators across India. Authorized 3B Scientific distributor in India. We provide complete support including pan-India delivery, on-site installation, faculty training workshops, and prompt after-sales maintenance. All our products carry genuine manufacturer warranty and meet ISO 9001 quality standards. **Frequently Asked Questions** Q: Is the Intradermal Injection Simulator compliant with NMC simulation lab requirements? A: Yes, this simulator supports the skill-based training competencies specified under NMC's CBME framework for Indian medical and nursing institutions. Q: Does SEM Trainers provide installation and training support? A: Absolutely. We offer on-site installation, hands-on demo sessions, and faculty orientation workshops for every product delivered, anywhere in India. Q: What is the warranty and how do I get support? A: All products come with standard manufacturer warranty. For support, reach us at [sem@semtrainers.com](mailto:sem@semtrainers.com) or call +91-88495 63724.