



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

I.V. & I.M. Injection Arm

Item No. LF00698U
Weight 9kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

- Realistic Flashback

- Replaceable Skins
- Intramuscular Injection Site complete with Fingerprints

Simple and easy to use, Life/form® Venipuncture and Injection Training Arms are the ideal way to teach venipuncture and injection techniques including starting IVs and introducing Over the Needle IV catheters. This simulator ends your problem of finding willing volunteers, and helps to avoid patient injury and liability while the students are learning new procedures. Realism begins with fine detail - right down to the fingerprints. But most important is the realistic feeling of puncture through the vinyl skin and latex veins. The skin actually rolls as you palpate the vein and the characteristic “pop” can be felt as the needle penetrates the vein. Veins are accessible at the antecubital fossa, along the forearm, and at the back of the hand, making it possible to practice venipuncture at any of the common sites. Practice of intramuscular injection is equally realistic. A foam pad simulates muscle tissue, and the boundaries are identified by natural bonelike landmarks. Life/form® Venipuncture and Injection Arms are designed with replaceable skins and veins. Complete replacement skin kits are available and easy to use. To extend the life of each skin/vein set, an aerosol sealant is available to seal punctures in the tubing and prevent leakage. Includes a Life/form® arm with skin, 22 gauge needle, butterfly set, infusion set, 3 cc syringe, 12 cc syringe, one pint Life/form® blood, fluid supply bags, instruction manual, and hard carry case. Fluid supply stand not included. Size: 35" x 12" x 6". CAUTION: Use only Life/form® Replacements. Do not use Conventional Latex Tubing. This will cause malfunction of the simulator.

Why Choose the I.V. & I.M. Injection Arm for Your Nursing Skills Lab in India?

The Life/form® I.V. & I.M. Injection Arm (Item No. LF00698U, 9 kg, Make: Nasco Healthcare, USA) is a comprehensive venipuncture and injection trainer, covering IV starts, Over-the-Needle IV catheter introduction, and intramuscular injection technique on a single realistic arm — solving the perennial problem of finding willing volunteers for student practice while eliminating patient injury and liability risk.

Realism runs down to the fingerprints molded into the intramuscular injection site. The vinyl skin rolls realistically as students palpate the vein, and a distinct "pop" is felt as the needle penetrates — genuine tactile feedback, not a stylized approximation. Veins are accessible at the antecubital fossa, along the forearm, and at the back of the hand, covering all common venipuncture sites. For intramuscular practice, a foam pad simulates muscle tissue, with natural bonelike landmarks defining the correct injection boundaries.

Skins and veins are replaceable, and an aerosol sealant is available to extend the life of each set by sealing punctures in the tubing.

Important: Use only genuine Life/form® replacement parts. Conventional latex tubing will cause the simulator to malfunction.

Includes: Life/form® arm with skin, 22 gauge needle, butterfly set, infusion set, 3cc syringe, 12cc syringe, one pint Life/form® blood, fluid supply bags, instruction manual, hard carry case.

Please note: a fluid supply stand is **not included** — order separately (see our Fluid Supply Stand, LF01022U).

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Can I use standard latex tubing as a replacement? A: No — only genuine Life/form® replacement parts should be used, or the simulator will malfunction.

Q: Is a fluid supply stand included? A: No — it's sold separately (LF01022U).

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.