



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



I.M. Injection Simulator

Item No. 1008941 [W30502/1]
Weight 0.45kg
MPN: M118-1

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

The student can use this injection simulator to practice both intramuscular and subcutaneous injections in the upper arm. The simulator can be attached easily using the stick-on tape. The student can find the injection site by feeling for the acromion. LED lights will show whether the injection has been inserted in the correct place to the right depth. The warning signal will sound if the injection is carried out incorrectly. Water only should be used for injecting. The delivery contents include the liquid pouch and the stick-on tape.

Why Choose the Wearable I.M./Subcutaneous Injection Trainer for Your Nursing Skills Lab in India?

This compact, wearable injection simulator (Item No. 1008941 [W30502/1], 0.45 kg, MPN: M118-1) lets students practice both intramuscular and subcutaneous injection technique on the upper arm, worn directly by a volunteer or instructor using simple stick-on tape — no full manikin required.

Students locate the correct injection site by palpating the acromion, the same landmark identification skill used in real clinical practice. LED lights confirm whether the injection was inserted at the correct site and to the correct depth, and an audible warning signal sounds if the injection is performed incorrectly — giving students immediate, objective feedback. **Water only should be used for injecting** with this device.

Its compact, wearable design makes it a practical option for skills stations, group practice, or portable training scenarios where a full-size manikin isn't practical.

Includes: Liquid pouch, stick-on tape.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of leading clinical skills training equipment in India, with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What fluid should be used with this device? A: Water only.

Q: How does it confirm correct technique? A: LED lights confirm correct site and depth; an audible alarm sounds if the injection is performed incorrectly.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495

63724 for pricing and availability.