



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Life/form® GERi™ / KERi™ Manikin IV Arm - Light

Item No. LF04080
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

Features replaceable skins and latex veins that actually roll when palpated. A realistic flashback confirms proper needle placement, and the shoulder can be used for intramuscular injection training. For use with all **GERi™** and **KERi™** manikins.

Why Choose the Life/form® GERi™ / KERi™ Manikin IV Arm — Light for Your Nursing Skills Lab in India?

The Life/form® GERi™ / KERi™ Manikin IV Arm – Light (Item No. LF04080, Make: Nasco Healthcare, USA) is a replacement IV arm designed specifically for use with GERi™ and KERi™ manikins, giving your existing manikin realistic venipuncture and IV training capability.

Compatibility note: This arm is designed for use with all GERi™ and KERi™ manikins. Please confirm your manikin model before ordering.

The arm features replaceable skins and latex veins that roll naturally when palpated, closely mimicking the feel of a real patient's vein. A realistic flashback confirms correct needle placement, giving students clear feedback during IV insertion practice. The shoulder can also be used for intramuscular injection training, extending the arm's use beyond venipuncture alone.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which manikins is this arm compatible with? A: It's designed for use with all GERi™ and KERi™ manikins.

Q: Does this arm support intramuscular injection practice as well as IV? A: Yes — in addition to venipuncture, the shoulder can be used for intramuscular injection training.

Q: How does the arm confirm correct needle placement? A: A realistic flashback occurs when the needle correctly enters the vein.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.