



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Gaumard® Advanced Multipurpose Venous Training Left Arm - Light

Item No. SB32870AL
Weight
Make: Gaumard Scientific, USA

[Read More](#)

SKU:SB32870AL

Categories: Injections and Punctures

Product Description

IV, IM, and sub-Q training Arm rotates 180° at deltoid region Subtle venous network in arm and hand Cephalic, basilic, antecubital, radial, and ulnar veins Realistic "pop" as needle enters vein Intramuscular injection site in deltoid area Subcutaneous injection areas on the volar side of the forearm and the lateral side of the upper arm Squeeze bulb to increase or decrease venous pressure Veins stand out or collapse Easily replaceable skin and veins Administration of medication by intravenous bolus Simulation of infusion technique Blood collection exercises with simulated blood Simulation of clenched fist and tourniquet position Includes simulated blood concentrate, pressure bulb, blood dispensing bag, spare arm skin, funnel, powder, instruction manual, smoked Lucite® base with stand, and soft nylon carry bag

Why Choose the Gaumard® Advanced Multipurpose Venous Training Left Arm — Light — for Your Nursing Skills Lab in India?

The Gaumard® Advanced Multipurpose Venous Training Left Arm, Light (Item No. SB32870AL, Make: Gaumard Scientific, USA) is a comprehensive arm trainer for IV, intramuscular, and subcutaneous injection technique. A subtle venous network across the cephalic, basilic, antecubital, radial, and ulnar veins gives students realistic site identification practice, and the arm rotates 180° at the deltoid for natural positioning during training.

A squeeze bulb lets instructors increase or decrease venous pressure on demand, making veins stand out or collapse — useful for simulating different patient presentations within the same session. Students feel a realistic "pop" as the needle enters the vein, and can practice the complete IV workflow: bolus medication administration, infusion technique, blood collection exercises with simulated blood, and proper clenched-fist and tourniquet positioning before venipuncture.

Beyond venous access, the arm includes a dedicated intramuscular injection site in the deltoid area and subcutaneous injection zones on the volar forearm and lateral upper arm — covering the three most common injection routes on one trainer. Skin and veins are easily replaceable as they wear from repeated use.

Includes: Simulated blood concentrate, pressure bulb, blood dispensing bag, spare arm skin, funnel, powder, instruction manual, smoked Lucite® base with stand, soft nylon carry bag.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Gaumard Scientific products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which veins are represented on this trainer? A: The cephalic, basilic, antecubital, radial, and ulnar veins.

Q: Can venous pressure be adjusted during training? A: Yes — a squeeze bulb lets instructors make veins stand out or collapse to simulate different patient presentations.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.