




SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Paddy - Small Injection Pad

Item No. SB51839
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:SB51839

Categories: Injections and Punctures

Product Description

Use to simulate intradermal, intramuscular, and subcutaneous injections. Removable epidermis layer on pad surface allows for the simulation of a "wheal" during intradermal injections. Skin puncture resistance and the resistance throughout the injection closely mimic the feel of a real injection. Drainage holes on pad bottom for use with fluid-filled syringes 2-1/2 in. L x 2-1/2 in. W x 1-3/4 in. H. Clean product this time — brand is clearly stated (Nasco Healthcare, USA), item number and dimensions are specific, and the description is accurate and complete. Only a minor field gap: **Weight is blank** on the live page (empty row, not even a placeholder number) — worth filling in from Nasco's spec sheet when convenient, but not urgent since it's not misleading like the 0.001 kg cases.

Paddy — Small Injection Pad for Intradermal, Intramuscular & Subcutaneous Practice

Paddy — Small Injection Pad (Item No. SB51839, 2-1/2 in. L × 2-1/2 in. W × 1-3/4 in. H) by Nasco Healthcare, USA is a compact, standalone injection training pad built to simulate three distinct injection techniques — intradermal, intramuscular, and subcutaneous — on a single small-footprint accessory.

Realistic technique feedback in a small package. The pad's removable epidermis layer allows students to see and feel an actual "wheal" form during intradermal injection practice, a detail that's easy to miss on larger, less detailed trainers. Skin puncture resistance and the resistance felt throughout the injection are designed to closely mimic a real injection, so students build accurate muscle memory before working with patients. Drainage holes on the pad's underside support use with fluid-filled syringes for a more realistic practice session.

Why it's useful for nursing and medical skills labs: Paddy's small size makes it practical for individual student stations or skills kits, rather than requiring a full-size manikin or arm just to practice injection technique — useful for high-volume basic skills training where every student needs repeated hands-on practice.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills trainers to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Which injection techniques can be practiced on Paddy? A: Intradermal, intramuscular, and subcutaneous injections — all on the same pad.

Q: Can this be used with fluid-filled syringes? A: Yes — drainage holes on the pad's underside support practice with fluid-filled syringes.

Q: Is the epidermis layer replaceable? A: The pad's removable epidermis layer allows for realistic "wheal" formation during intradermal practice. Confirm with our team on replacement part availability if you're planning high-volume use.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.