



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Patrick - Large Injection Pad

Item No. SB51840
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:SB51840

Categories: Injections and Punctures

Product Description

- Use to simulate intramuscular, intradermal (when used with removable epidermis layers), subcutaneous, and Z-track injections
- Realistic puncture resistance mimics the feel of an actual injection
- Bottom drainage allows the pad to be used with fluid-filled syringes
- Pad skin can be grasped for more realistic injection scenarios
- 6 in. L x 6 in. W x 1-3/4 in. H

Patrick — Large Injection Pad for Multi-Technique Injection Training

Patrick — Large Injection Pad (Item No. SB51840, 6 in. L × 6 in. W × 1-3/4 in. H) by Nasco Healthcare,

USA is a standalone injection training pad built to support four distinct injection techniques: intramuscular, subcutaneous, Z-track, and intradermal (when used with removable epidermis layers).

A wider technique range than a basic pad. Where a smaller pad typically covers the core three injection types, Patrick adds Z-track injection — a technique used for medications that can irritate or stain skin if they leak back along the needle path — giving students a more complete injection-technique practice surface. Realistic puncture resistance mimics the feel of an actual injection, and the pad's larger size means the skin can be grasped during practice, replicating how a nurse or clinician would stabilize tissue on a real patient for more realistic scenario training.

Practical for repeated use: bottom drainage allows the pad to be used with fluid-filled syringes for added realism, and the larger 6 × 6 inch surface gives more room for repeated needle sticks across a training session before the pad needs replacing.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorised distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills trainers to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What injection techniques can be practiced on Patrick? A: Intramuscular, subcutaneous, Z-track, and intradermal injection (intradermal requires removable epidermis layers — confirm with our team whether these are included or need to be ordered separately).

Q: How is this different from the smaller Paddy pad? A: Patrick is larger (6 × 6 in. vs. Paddy's 2-1/2 × 2-1/2 in.) and adds Z-track injection technique, plus a grabbable pad skin for more realistic scenario practice.

Q: Can this be used with fluid-filled syringes? A: Yes — bottom drainage supports fluid-filled syringe use.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.