



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Multi Energy CT Quality Assurance Phantom (TR-I)

Item No. PH-75B
Weight 23 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Water Equivalent Material, various inserts and empty bottles enable to verify the appropriate Multi-Energy CT settings

Features

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Optional parts

Related products

Production ? Development Supervision

Update

1. Phantom using innovative Water Equivalent Material /
2. Empty bottles enable to put various items and check how they react to Multi Energy CT / 3. Save time and efforts to produce custom-made water phantoms / 4. Two different sizes of body. (TR-I, TR-J)

ME-CT image analysis protocol / Metal artifact reduction / Reduction of contrast media

1 truck phantom / 1 internal cylindrical phantom / 9 filling rods for dosimeter holes / 8 WEM inserts / 3 iodine inserts (4, 8, 12mgI/mL) / 1 titanium insert / 1 soft tissue insert / 1 water container/ 20 empty sample bottles/ 8 spacers for sample bottles/ manual

□36.3×D18×H26.2cm

W36.3 x D26.2 x H18.0 cm/W14.3 x D10.3 x H7.1 in

23kg

W44 x D36 x H42cm / W17.3 x D14.2 x H16.5in

Polyurethane

Angle adjustment holder

[Multi Energy CT Quality Assurance Phantom \(TR-J\)](#)

Katsuhiro Ichikawa, Ph.D., / Institute of Medical, / Pharmaceutical and Health Sciences, / Kanazawa University

October 25, 2024

Why Choose the Multi Energy CT Quality Assurance Phantom (TR-I) for ME-CT Protocol Evaluation in India?

The Multi Energy CT Quality Assurance Phantom TR-I (Item No. PH-75B, Brand: Kyoto Kagaku, Japan) uses innovative Water Equivalent Material construction, with empty sample bottles that let you insert and test how various materials respond under multi-energy CT settings — saving the time and effort of producing custom-made water phantoms. This phantom is available in two body sizes (TR-I and TR-J). Developed under the supervision of Katsuhiro Ichikawa, Ph.D., Institute of Medical, Pharmaceutical and Health Sciences, Kanazawa University.

Training applications:

- ME-CT (multi-energy CT) image analysis protocol development
- Metal artifact reduction evaluation
- Contrast media reduction assessment

Set includes: 1 truck phantom, 1 internal cylindrical phantom, 9 filling rods for dosimeter holes, 8 WEM inserts, 3 iodine inserts (4, 8, 12 mgI/mL), 1 titanium insert, 1 soft tissue insert, 1 water container, 20 empty sample bottles, 8 spacers for sample bottles, manual.

Specifications:

- Size approx. W36.3×D18×H26.2cm
- Packing size approx. W44×D36×H42cm (W17.3×D14.2×H16.5in)
- Weight approx. 23kg
- Materials: Polyurethane

Optional accessory: Angle adjustment holder.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying multi-energy CT quality assurance equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What can the empty sample bottles be used for? A: They let you insert your own test materials to see how they respond under multi-energy CT settings — useful for custom research without needing to build a separate water phantom.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.