



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com







DVT (Deep Vein Thrombosis) Phantom

Item No. US-21
Weight 3 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

For ultrasound screening of DVT (Deep Vein Thrombosis)

Features

Training skills / Applications

Set includes

Size (approx.)

Weight (approx.)

Production ? Development Supervision

Update

Faithfull reproduction of blood vessel meandering, branching, and thrombus. The built-in blood vessels allows for milking method as well as observation of compression of the vessels by applying pressure with a transducer. Familiarize with assessment of venous thrombosis in the lower limbs.

Screening for deep venous thrombosis in the lower limbs
Observation of blood flow by milking method

1 lower limb manikin 1 ultrasound calf phantom (incorporated) 1 storage case 1 instruction manual

W27 x D17 x H47 cm

3.0kg

Kenji Karino, Professor, Clinical Skill-Up Center, Simane University Hospital

May 16, 2023

Why Choose the DVT (Deep Vein Thrombosis) Phantom for Your Vascular Ultrasound Training Program in India?

The DVT (Deep Vein Thrombosis) Phantom (Item No. US-21, 3.0kg, Brand: Kyoto Kagaku, Japan) is purpose-built for ultrasound screening training of deep vein thrombosis. Production and development were supervised by Kenji Karino, Professor at the Clinical Skill-Up Center, Shimane University Hospital.

The phantom faithfully reproduces blood vessel meandering, branching, and thrombus formation, with built-in blood vessels that allow both the milking method and observation of vessel compression by applying pressure with a transducer — the two core techniques used in real clinical DVT screening. This lets trainees build genuine competency in assessing venous thrombosis in the lower limbs before working with real patients.

Training applications: Screening for deep venous thrombosis in the lower limbs, observation of blood flow by milking method.

Set includes: Lower limb manikin, ultrasound calf phantom (incorporated), storage case, instruction manual.

Specifications: Size W27×D17×H47cm, weight 3.0kg.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and vascular ultrasound training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What techniques can be practiced? A: The milking method and observation of vessel compression by transducer pressure — the two core clinical DVT screening techniques.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.