



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Intravesical Urine Volume Measurement Phantom

Item No. US-16
Weight 11.7 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Pocket-sized ultrasound is gaining popularity and starting to be utilized more in daily care for examination. This method is preferred as catheterization may cause urinary infection.

Features

4 variations of interchangeable inserts feature different patient scenarios.

Training skills / Applications

Case / Pathology

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Related products

Production ? Development Supervision

CE certified/CE Certification?

Update

Handling and manipulation of transducers on bladder ultrasound/ Scanning the bladder to measure urine volume and findings.
Urine volume 50mL, 150mL, 300mL, Urinary retention + balloon catheter
1 lower torso manikin, 4 ultrasound blocks, 1 talcum powder, 1 strage case, 1 instruction manual
30 x 26 x 18 cm 12 x 10 x 7 in
W52 x D49 x H29cm / W20.5 x D19 x H11.5in
Phantom: 1.5kg / 3.3lbs / Body: 3.2kg / 7.1lbs
11.7kg/25.8lb
Soft resin, hard resin, Latex free
MW2A Male Catheterization/Enema Simulator MW52 Female Catheterization/Enema Simulator US-20 Ultrasound Guided Tube Feeding Phantom
Dr. Tadashi Kobayashi, Assistant Professor of Hirosaki University Hospital
No
November 1, 2020