



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




NEW! Ultrasound Guided Tube Feeding Phantom

Item No. US-20
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Anatomically correct neck model for Ultrasound Guided Tube Feeding Training

Features

Training skills / Applications

Set includes

Size (approx.)

Weight (approx.)

Related products

CE certified/CE Certification?

Update

1. Anatomically correct neck model Thyroid Gland Trachea Esophagus Sternum Collarbone Carotid Artery Jugular Vein 2. Realistic neck ultrasound images 3. Usable with real ultrasound diagnostic devices 4. Allows the confirmation of the proper Nasogastric Tube Placement with ultrasound 5. Inner switch for tube insertion direction: trachea/esophagus for tube location image recognition tests (recognition by ultrasound of proper esophagus insertion or wrong tracheal insertion)

| Proper handling and manipulation of probe on neck ultrasound | Feeding tube insertion | Confirmation of tube placement by ultrasound | Neck scan to evaluate the tube position.

Main Body Ultrasound neck block Nasogastric Tube Syringe Talcum powder Storage case

W29 x D18 x H17cm(Body) W17 x D10 x H10cm(cube) W11.4 x H6.7 x D7.1inch

1.8kg(body) / 1.1kg(cube)

[MW8 Tube Feeding Simulator\(NG, OG, and PEG\)](#)

No

November 1, 2020

Why Choose the Ultrasound Guided Tube Feeding Phantom for Your Point-of-Care Ultrasound Training in India?

The Ultrasound Guided Tube Feeding Phantom (Item No. US-20, Brand: Kyoto Kagaku, Japan) is an anatomically correct neck model purpose-built for training ultrasound-guided nasogastric tube placement confirmation — a faster, radiation-free alternative to X-ray verification.

The model reproduces the thyroid gland, trachea, esophagus, sternum, collarbone, carotid artery, and jugular vein, and is usable with real ultrasound diagnostic devices for genuinely realistic neck ultrasound images. An internal switch lets instructors set the tube's simulated path to either the trachea or esophagus, so students can practice the critical clinical skill of recognizing correct esophageal tube placement versus incorrect tracheal insertion purely through ultrasound image interpretation — the judgment that matters most in real patient care.

Training applications: Proper handling and manipulation of the probe on neck ultrasound, feeding tube insertion, confirmation of tube placement by ultrasound, and neck scanning to evaluate tube position.

Set includes: Main body, ultrasound neck block, nasogastric tube, syringe, talcum powder, storage case.

Specifications: Body size W29×D18×H17cm, cube size W17×D10×H10cm (W11.4×H6.7×D7.1in). Body weight 1.8kg, cube weight 1.1kg.

Please note: this device is **not CE certified** — confirm with your institution whether this affects your procurement requirements.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and point-of-care ultrasound training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: Can this be used with real ultrasound equipment? A: Yes — the model is designed for use with real ultrasound diagnostic devices.

Q: What does the internal switch feature do? A: It sets the tube's simulated path to trachea or esophagus, letting students practice recognizing correct versus incorrect placement by ultrasound alone.

Q: Is this CE certified? A: No — please confirm with your institution whether this affects your procurement requirements.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.