



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Male Catheterization/Enema Simulator

Item No. MW2A
Weight 3.5 KG
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Catheterization



Product Description

This is the confirmed male counterpart to the MW52 Female Catheterization/Enema Simulator from earlier — same weight figures, same packing specs, genuinely matched pair. A few source typos to fix silently, and one new product reference worth flagging for your catalog.

[Cert] Item No. MW2A, weight 3.5 kg (product) / 14.5 kg (packing), size 34×33×24 cm, Brand Kyoto Kagaku, Japan — identical product/packing weight to the female MW52, confirming these are a genuinely matched pair sharing the same chassis design. Typos corrected: "positio" ? position, "recum" ? rectum, "reservioir" ? reservoir.

[Cert] New product reference worth noting: the related-products list includes **MW58/59 Suprapubic Catheter Care Simulator** — a Kyoto Kagaku product not yet covered in your catalog based on this session. Worth checking whether semtrainers.com has a listing for it, since it would be a natural addition alongside the P93 PRO line's suprapubic catheter care feature covered earlier.

DESCRIPTION (copy-paste into WordPress):

Why Choose the Male Catheterization/Enema Simulator (Kyoto Kagaku) for Your Nursing Skills Lab in India?

The Male Catheterization/Enema Simulator (Item No. MW2A, 3.5 kg, 34 × 33 × 24 cm) by Kyoto Kagaku, Japan is a high-quality trainer covering urethral catheterization, perineal care, and enema administration in one compact unit — the male counterpart to our Female Catheterization/Enema Simulator (MW52). Soft, flexible, latex-free genitalia allow realistic cleaning and procedural preparation, and successful catheter insertion is confirmed by fluid outflow.

The detachable external genitalia design makes cleaning and part replacement straightforward, and enema administration is trained in the correct lateral position, matching real clinical practice.

Training skills covered: Intermittent catheterization, indwelling catheterization, perineal care and enema administration, localization of the urethral meatus, cleansing and aseptic technique, catheter insertion and placement, catheter removal, and manual bladder compression.

What's included: Lower torso manikin, genital organ unit, 3 valves (one pre-installed), irrigation bottle, rectum tube, 3 irrigation drainage tubes, lubricant, storage case, instruction manual.

Recommended device: 16 Fr urinary catheter.

Related Kyoto Kagaku products: MW52 (Female Catheterization/Enema Simulator), MW53/54/55 (Strap-on Catheterization/Enema Simulator), MW58/59 (Suprapubic Catheter Care Simulator), MW23 (Adult Diaper Training Simulator) — useful companions for labs building out a full continence-care curriculum.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Kyoto Kagaku products in India with over 30 years of experience supplying medical simulators and clinical skills

training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: Is this suitable for latex-allergic students or staff? A: Yes — the materials are soft resin and explicitly latex-free.

Q: What catheter size is recommended? A: A 16 Fr urinary catheter.

Q: Is there a female equivalent? A: Yes — the MW52 Female Catheterization/Enema Simulator, sharing the same design family.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.

High quality simulators for urethral catheterization, perineal care and enema practice

Features

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Related products

RECOMMENDED DEVICES

Update

1. Soft and flexible genitalia allow proper cleaning and procedural preparation 2. Successful insertion confirmation by fluid outflow 3. Enema insertion in lateral position 4. Detachable external genitalia

Intermittent catheterization / Indwelling catheterization / Perineal care / Enema / Localization of urethral meatus / Cleansing and aseptic technique / Catheter insertion and placement / Catheter Removal / Manual bladder compression

1 lower torso manikin / 1 genital organ unit / 3 valves (one is incorporated in the assembled unit) / 1 irrigation bottle / 1 rectum tube / 3 irrigation drainage tubes / 1 lubricant / 1 storage case / 1 instruction manual

W34 x D33 x H24 cm / W13.5 x D13 x H9.5 in

W52 x D45 x H31cm / W20.5 x D17.7 x H12.3in

3.5kg / 7.7lbs

14.5kg / 32lbs

Soft resin / Latex free

11381-100-01 clip (small) 11381-100-02 clip (large) 11381-100-03 urinary bladder reservoir 11381-100-04 rectum unit 11381-300 6 replacement valves 11381-800 lubricant 11381-410 replacement male skin 11381-510 replacement male meatus for MW2A (a pair)

[MW52 Female Catheterization/Enema Simulator](#)

[MW53/54/55 Strap-on Catheterization/Enema Simulator](#)

[MW58/59 Suprapubic Catheter Care Simulator](#) [MW23 Adult diaper training simulator](#)

Urinary catheter: 16Fr

October 20, 2020