



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Smart Stat Xtreme Trauma Bleeding Arm

Item No. 101-481
Weight 11 lbs
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Moulage and Wound Simulation

Product Description

These graphic explosion wounds offer detailed enhancements to challenge any practitioner's level of care. These moulages feature missing muscle and bone, flesh avulsions, fatty tissue, and bone ends. The blood supply line from the SMART STAT simulator connects quickly, without blood splatter, to the moulage, and has pressurized blood flow to create a messy wound. Proper application of a tourniquet will halt the blood flow. Xtreme Bleeding Moulages can be used with the STAT line of patient simulators. External blood supplies must be used if not attached to the SMART STAT. Xtreme Trauma Bleeding Arm. This arm can replace the left arm of the SMART STAT and consists of a molded moulage piece with a structure that offers a simple connection to the existing IM injection site of a SMART STAT patient simulator. Shipping weight 11 lbs.

Smart Stat Xtreme Trauma Bleeding Arm (Nasco Healthcare) — Pressurized

Bleeding Moulage for STAT Simulators in India

The Smart Stat Xtreme Trauma Bleeding Arm (Item No. 101-481, Brand: Nasco Healthcare, USA) is a graphic explosion-wound moulage designed to challenge advanced levels of trauma care, featuring missing muscle and bone, flesh avulsions, fatty tissue, and visible bone ends.

The blood supply line from the SMART STAT simulator connects quickly to the moulage without blood splatter, delivering pressurized blood flow to create a realistic, messy wound. A properly applied tourniquet will halt the blood flow, giving trainees a genuine feedback loop for hemorrhage control technique. This arm replaces the left arm of the SMART STAT simulator and connects to the existing IM injection site via a simple molded moulage structure. An external blood supply is required if the moulage is not attached directly to a SMART STAT unit.

Specifications: Weight (ship) 11 lbs (approx. 5 kg). Brand: Nasco Healthcare, USA. Compatible with the STAT line of patient simulators.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying advanced trauma and hemorrhage-control training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

FAQ Q: How does the pressurized bleeding feature work? **A:** The SMART STAT blood supply line connects to the moulage without splatter, delivering pressurized flow that a correctly applied tourniquet will stop — giving trainees real feedback on hemorrhage control technique.

Q: How does this arm attach to the simulator? **A:** It replaces the left arm of a SMART STAT patient simulator, connecting via the existing IM injection site structure.

Q: How do I order from SEM Trainers? **A:** Contact sem@semtrainers.com or call +91-88495 63724 for pricing and availability.